

Skin Markers

Pre-Op Mini Marker

Features

- Assists in the prevention of wrong site/wrong surgery
- Less waste than full size skin marker
- Eliminates opportunity for cross-contamination by reducing the temptation to reuse
- Comfortable thickness and length for ease of use
- Convenient storage - 100 each tub
- FDA approved ink for skin - non-tattooing
- Less expensive than a standard skin marker
- Each marker is labelled 'single patient use only'
- TGA listed - ARTG#154608



Product

Pre-Op Mini Markers

Pre-Op Mini Markers

Catalogue No. Packaging

DRY26-017NS 800/Carton

DRY26-018NS 200/Carton

Sterile Skin Markers

Features

- For smooth, clean line markings on patient's skin
- Prep-resistant genetian violet ink
- Non-toxic and non-tattooing
- Each barrel legibly ruled in centimeters
- Flexible 15cm ruler allows measurement of body curves
- Individually packed in sterile peel pouches for fibre-free introduction into sterile field
- Dispenser boxes minimise storage and maximise convenience
- Latex Safe



Product

Sterile Skin Markers, standard tip w/ ruler; blue barrel

Sterile Skin Markers, fine tip w/ ruler; white barrel

Catalogue No. Packaging

DRY26-001 25/Carton

DRY26-002 25/Carton