




# Solutions for Oxygen Therapy & Airway Management

Designed to ensure optimal performance,  
safety and patient comfort

**ThermoFisher**  
SCIENTIFIC



Thermo Fisher Scientific is pleased to provide Healthcare Professionals with our comprehensive range of oxygen therapy and airway management products.

Our promed® respiratory support and airway management range has been designed to ensure optimal performance and safety while ensuring patient comfort during therapy.

Our comprehensive range of respiratory support products allows for the provision of supplemental oxygen, while the airway management range ensures safe and effective management of airways and containment of fluids.

All of the respiratory support and airway management products offered have been validated for use in the Australian healthcare environment and certified by the Therapeutic Goods Administration (TGA).

We pride ourselves on the design and quality of our products that have become a promed® hallmark for more than 40 years; we will continue to ensure we exceed exceptional levels of safety, functionality and reliability with competitive pricing.

# Contents



Oxygen Masks	4
--------------	---

Aerosol Masks	5
---------------	---

Venturi Mask Kit	5
------------------	---

Anaesthetic Masks	5
-------------------	---

Oxygen Tubing	6
---------------	---



Aerosol Tubing — Corrugated	6
-----------------------------	---

Nasal Cannulae	7
----------------	---

Nebulisers & Kits	8
-------------------	---

Nasopharyngeal Airways	8
------------------------	---

Oropharyngeal Airways	9
-----------------------	---

Laryngeal Mask Airways	10
------------------------	----



Intubating Stylets & Oxygen Catheters	11
--	----

Respiratory Adaptors & Anaesthesia Accessories	12
---	----

Incentive Spirometer	13
----------------------	----

Yankauer Suction Devices	13
--------------------------	----

Suction Tubing	13
----------------	----

# promed Oxygen Masks

The promed Oxygen Masks are made from high quality materials which result in a soft, clear finish designed to optimise comfort and facilitate easy monitoring of patient respiratory status.

Adult and paediatric masks feature an adjustable nose clip and all masks come with adjustable elastic braid strip to ensure a customised fit for every patient. All promed masks are compatible with standard oxygen tubing connector.

## promed Medium Concentration Masks

Item No.	Product	Qty
PMD06006	Med Conc Mask Adult Standard No Tube	50/pack
PMD06005-GC	Med Conc Mask Adult Standard 3m Oxygen Tubing	50/pack
PMD06000-GC	Med Conc Mask Adult Elongated No Tube	50/pack
PMD06104-GC	Med Conc Mask Adult Elongated 3m Tubing No-Crush Tubing	50/pack
PMD06007-GC	Med Conc Mask Paediatric Elongated No Tube	50/pack
PMD06105-GC	Med Conc Mask Paediatric Elongated 2.1m Tubing No-Crush Tubing	50/pack
PMD06010	Med Conc Mask Infant O <sub>2</sub> Tube Connector 2.1m O <sub>2</sub> Tubing	50/pack

## promed High Concentration Masks

Item No.	Product	Qty
PMD06001	Hi Conc Mask Rebreather with Tubing 1 Valve Elongated	50/pack
PMD06009-GC	Hi Conc Mask Paediatric O <sub>2</sub> Tube 1 Valve Non-Rebreather 3m	50/pack
PMD06002	Hi Conc Mask Non Rebreather with Tubing 1 Valve Elongated	50/pack
PMD06003-GC	Hi Conc Mask Non Rebreather Tube 2 Valves Elongated	50/pack
PMD06008	Hi Conc Mask Paediatric Rebreather No Tube	50/pack
PMD06106	Hi Conc Mask Paediatric Non Rebreather Tube 2 Valves Elongated	50/pack

## promed Tracheostomy Masks

Item No.	Product	Qty
PMD06011	Tracheostomy Mask, Adult No Tubing	50/pack
PMD06012	Tracheostomy Mask, Paediatric No Tubing	50/pack

## promed 3-in-1 Emergency Mask

Item No.	Product	Qty
PMD06004	Emergency Mask 3-In-1 Elongated with 2.1m Oxygen Tubing	50/pack



## promed Aerosol Masks

The promed Aerosol Masks are designed to deliver aerosol therapy while optimising patient comfort. Made from high quality materials, adult and paediatric masks feature adjustable nose clips for custom fitting and an adjustable elastic braid strap. The promed® masks connector accepts standard 22mm aerosol tubing.

Item No.	Product	Qty
PMD02000	Aerosol Masks Adult, Elongated (Under Chin) No Tubing	50/pack
PMD02001	Aerosol Masks Adult, Face Tent, No Tubing	50/pack
PMD02003	Aerosol Masks Adult, Standard (Under Lip), No Tubing	50/pack
PMD02004	Aerosol Masks Paediatric, Elongated (Under Chin), No Tubing	50/pack



## promed Venturi Mask Kit

The promed Venturi Mask kits allow for delivery of fixed concentration of oxygen. The fixed delivery of oxygen is achieved through the use of either colour-coded diluters or an adjustable diluter. Both kits are conveniently pre-packaged.

Item No.	Product	Qty
PMD02002	Venturi Mask Kit Adult with Adjustable Diluter No-Crush Lumen Oxygen Tubing	50/pack
PMD02005	Venturi Mask Kit Adult All Diluters No-Crush Lumen Oxygen Tubing	50/pack



## promed ClearVue Plus Anaesthetic Masks

promed ClearVue Plus face masks have a clear cushion for easy viewing and monitoring of patient condition. The cushion features a syringe inflation port at the bridge of the mask for greater inflation convenience. Masks are available in five sizes. The larger sizes come fitted with collar rings which can be removed for hand held procedures.

Item No.	Product	Qty
PMD03010	ClearVue Plus Mask Infant OD 15mm Size 1 No Ring	20/pack
PMD03020	ClearVue Plus Mask Child ID 22mm Size 2	20/pack
PMD03030	ClearVue Plus Mask Adult Small ID 22mm Size 3	20/pack
PMD03040-GC	ClearVue Plus Mask Adult Medium ID 22mm Size 5	20/pack
PMD03050-GC	ClearVue Plus Mask Adult Large ID 22mm Size 6	20/pack



## promed Oxygen Tubing

The promed Oxygen Tubing products are manufactured from flexible non-toxic material. They are available in both standard or Non-Crush Lumen. Standard Lumen has a smooth internal bore while the No-Crush Lumen has a star-shaped interior cross section which reduces the likelihood of interruption of oxygen flow should the tubing be compressed or kinked.

Item No.	Product	Qty
PMD06100-GC	Oxygen Tubing 3m Standard	50/pack
PMD06101	Oxygen Tubing 2.1m No-Crush Lumen	50/pack
PMD06111-GC	Oxygen Tubing 3m No-Crush Lumen	50/pack
PMD06102-GC	Oxygen Tubing 4.2m No-Crush Lumen	50/pack
PMD06103	Oxygen Tubing 7.6m No-Crush Lumen	50/pack
PMDI10732	Oxygen Bubble Tubing 30m	1/pack



## promed Aerosol Tubing — Corrugated

The promed aerosol corrugated tubing is flexible and lightweight for maximum performance. Available in a variety of lengths with cuttable cuffs every 15cm (6inch) for convenient length adjustments.

Item No.	Product	Qty
PMD02015	Aerosol Tubing 1.5cm (60inch) with Cuffs every 15cm (6inch)	50/pack
PMD02014	Aerosol Tubing 1.8cm (72inch) with Cuffs every 15cm (6inch)	50/pack
PMD02016	Aerosol Tubing 30.5m (100ft) with Cuffs every 15cm (6inch)	50/pack




# promed Nasal Cannulae

The promed Nasal Cannulae are available in two different styles: Soft-Touch and Standard. The soft-touch Nasal Cannula range is made using a special manufacturing process which results in a flexible soft finish for increased patient comfort. They are light weight, fit comfortably over the ear and are available with straight tip for infants and flared curved tip in adult and paediatric sizes. Standard Cannula are available in a variety of configurations which include oxygen tubing with or without No-Crush Lumen.

## promed Standard Nasal Cannula

Item No.	Product	Qty
PMD08001-GC	Adult, Straight Tip, No Tubing	50/pack
PMD08004	Adult, Flared Curved TIP, No Tubing	50/pack
PMD08000-GC	Adult, Straight Tip with 3m Tubing	50/pack
PMD08005-GC	Adult, Flared Curved Tip with 3m Tubing	50/pack
PMD08006-GC	Adult, Flared Curved Tip with 3m Tubing, No-Crush Lumen	50/pack
PMD08007	Adult, Flared Straight Tip with 2.1m Tubing No-Crush Lumen	50/pack



## promed Soft-Touch Nasal Cannula

Item No.	Product	Qty
PMD08003	Adult, Curved Tip, No Tubing	50/pack
PMD08002-GC	Adult, Curved Tip 3m Tubing	50/pack
PMD08008	Adult, Curved Tip 2.1m Tubing No-Crush Lumen	50/pack
PMD08009	Paediatric Curved Tip 2.1m No-Crush Lumen	50/pack
PMD08010	Infant Straight Tip. 2.1m No-Crush Lumen	50/pack



## promed Et CO<sub>2</sub> MAC-SAFE Nasal Cannula

Item No.	Product	Qty
PMD08011	Adult 2.1m Tubing Female Luer, Fit-All Connector, 10cm/4 inch Gas Sampling Line	20/pack
PMD08013	Adult 2.1m Tubing Female Luer, Fit-All Connector, Gas Sampling Line	20/pack
PMD08012-GC	Adult 3m Tubing Male Luer, Fit-All Connector, Gas Sampling Line	20/pack



## promed Nebulisers & Kits

The promed Econi-Mist II and New OPTI-MIST PLUS are small volume nebulisers with easy to use designs. They are spill-proof and provide fast and effective medication delivery. Medication capacity is 5cc/mL with measurements conveniently graduated on the nebuliser medication cups for ease of use. The Large Volume Nebuliser has a dial design measuring accurate O<sub>2</sub> levels from 35% to 100%. Available in a 400cc/mL bottle with minimum and maximum level markings clearly defined. Connector accepts standard 22mm aerosol tubing.



Item No.	Product	Qty
PMD02600	Opti-Mist Plus Nebuliser, 2.1m Tubing T-Piece	20/pack
PMD02601	Econo-Mist II Nebuliser System with Opti-Mist Plus Nebuliser	50/pack
PMD02602-GC	Opti-Mist Plus Nebuliser 2.1m Tubing Adult Aerosol Mask, Elongated (Under Chin)	50/Pack
PMD02603-GC	Opti-Mist Plus Nebuliser, 2.1m Tubing Paediatric Mask Elongated (Under Chin)	50/pack
PMD02604	Econ-Mist II Nebuliser System, with Opti-Mist Plus Nebuliser	20/pack
PMD02605-GC	Opti-Mist Nebuliser Only	50/pack
PMD02606	Opti-Mist Plus Nebuliser Kit, 2.1m Tubing T-Piece, Mouth Piece	50/pack
PMD02013	Large Volume Nebuliser, 400cc/ml with Adaptor	50/pack



## promed Nasopharyngeal Airways

A nasopharyngeal airway is a short tube which is designed to be inserted into the nasal passageway to secure an open airway. It can be used safely on unconscious patients, and can also provide access for suction, both in the field and hospital environments.

- Single use
- Soft and pliable airway with smooth and rounded bevelled tip for increased patient comfort safety
- Resists collapse or kinking to help ensure a patent airway
- Size ranging from 20 to 33 French Gauge (FR)
- Latex and DEHP Free



Item No.	Product	Qty
PMD03112	Nasopharyngeal Airway PVC, ID 4.0mm, OD 6.5mm (20FR)	20/pack
PMD03114	Nasopharyngeal Airway PVC, ID 4.5, OD 7.0mm (21FR)	20/pack
PMD03116	Nasopharyngeal Airway PVC, ID 5.5, OD 8.0mm (24FR)	20/pack
PMD03118	Nasopharyngeal Airway PVC, ID 6.0, OD 8.5mm (26FR)	20/pack
PMD03119	Nasopharyngeal Airway PVC, ID 6.5, OD 9.0mm (27FR)	20/pack
PMD03120	Nasopharyngeal Airway PVC, ID 7.0, OD 9.5mm (29FR)	20/pack
PMD03122	Nasopharyngeal Airway PVC, ID 7.5, OD 10.0mm (30FR)	20/pack
PMD03124	Nasopharyngeal Airway PVC, ID 8.0, OD 10.5mm (32FR)	20/pack
PMD03126	Nasopharyngeal Airway PVC, ID 5.5, OD 11.0mm (33FR)	20/pack

# promed Oropharyngeal Airways

The Guedel Airway is designed to maintain a secure Oropharyngeal airway for the unconscious patient in emergency pre-hospital care, within the hospital environment and during the anaesthetic peri-operative period.

- Rounded atraumatic edges
- Smooth airway path for easy airway suctioning
- Bite block to prevent biting of tongue or endotracheal tube
- Colour coded for easy identification of sizes
- Single packed/sterile



Item No.	Product	Qty
PMD04100	Guedel Airway Size 000 (40mm) Pink	100/pack
PMD04000	Guedel Airway Size 00 (50mm) Dark Blue	100/pack
PMD04001	Guedel Airway Size 0 (60mm) Black	100/pack
PMD04010	Guedel Airway Size 1(70mm) White	100/pack
PMD04020-GC	Guedel Airway Size 2 (80mm) Green	50/pack
PMD04030-GC	Guedel Airway Size 3 (90mm) Yellow	50/pack
PMD04040	Guedel Airway Size 4 (100mm) Red	50/pack
PMD04050	Guedel Airway Size 5 (110mm) Light Blue	50/pack
PMD04060	Guedel Airway Size 6 (120mm) Orange	50/pack

# promed Laryngeal Mask Airways

The promed Laryngeal Mask Airways (LMA) are available as single-use and reusable, both in medical-grade PVC and silicone. The reusable LMA can be utilised and sterilised up to 40 times. The LMAs offer a low pressure seal around the laryngeal mask, around the laryngeal inlet and permit gentle positive pressure and ventilation. The reinforced LMA is available with spiral stainless steel for reinforcement. All promed LMAs feature an X-ray detectable line on the tube for patient safety.



Item No.	Product	Qty
PMDL100PS	Single Use, Laryngeal Mask PVC Size 1	1/pack
PMDL150PS	Single Use, Laryngeal Mask PVC Size 1.5	1/pack
PMDL200PS	Single Use, Laryngeal Mask PVC Size 2.0	1/pack
PMDL250PS	Single Use, Laryngeal Mask PVC Size 2.5	1/pack
PMDL300PS	Single Use, Laryngeal Mask PVC Size 3	1/pack
PMDL400PS	Single Use, Laryngeal Mask PVC Size 4	1/pack
PMDL500PS	Single Use, Laryngeal Mask PVC Size 5	1/pack
PMDL100SS	Single Use, Laryngeal Mask Silicone Size 1	1/pack
PMDL150SS	Single Use, Laryngeal Mask Silicone Size 1.5	1/pack
PMDL200SS	Single Use, Laryngeal Mask Silicone Size 2	1/pack
PMDL250SS	Single Use, Laryngeal Mask Silicone Size 2.5	1/pack
PMDL300SS	Single Use, Laryngeal Mask Silicone Size 3	1/pack
PMDL400SS	Single Use, Laryngeal Mask Silicone Size 4	1/pack
PMDL500SS	Single Use, Laryngeal Mask Silicone Size 5	1/pack



## promed Laryngeal Mask Airways, Single Use, Reinforced

Item No.	Product	Qty
PMDL100SRS	Single Use, Reinforced Laryngeal Mask Silicone Size 1	1/pack
PMDL150SRS	Single Use, Reinforced Laryngeal Mask Silicone size 1.5	1/pack
PMDL200SRS	Single Use, Reinforced Laryngeal Mask Silicone size 2	1/pack
PMDL250SRS	Single Use, Reinforced Laryngeal Mask Silicone Size 2.5	1/pack
PMDL300SRS	Single Use, Reinforced Laryngeal Mask Silicone size 3	1/pack
PMDL400SRS	Single Use, Reinforced Laryngeal Mask Silicone size 4	1/pack
PMDL500SRS	Single Use, Reinforced Laryngeal Mask Silicone size 5	1/pack
PMDL100SRP	Single Use, Reinforced Laryngeal Mask PVC size 1	1/pack
PMDL150SRP	Single Use, Reinforced Laryngeal Mask PVC size 1.5	1/pack
PMDL200SRP	Single Use, Reinforced Laryngeal Mask PVC size 2	1/pack
PMDL250SRP	Single Use, Reinforced Laryngeal Mask PVC size 2.5	1/pack
PMDL300SRP	Single Use, Reinforced Laryngeal Mask PVC size 3	1/pack
PMDL400SRP	Single Use, Reinforced Laryngeal Mask PVC size 4	1/pack
PMDL500SRP	Single Use, Reinforced Laryngeal Mask PVC size 5	1/pack



### promed Laryngeal Mask Airways Reusable Standard

Item No.	Product	Qty
PMDL100SR	Reusable Laryngeal Mask Silicone Size 1	1/pack
PMDL150SR	Reusable Laryngeal Mask Silicone Size 1.5	1/pack
PMDL200SR	Reusable Laryngeal Mask Silicone Size 2	1/pack
PMDL250SR	Reusable Laryngeal Mask Silicone Size 2.5	1/pack
PMDL300SR	Reusable Laryngeal Mask Silicone Size 3	1/pack
PMDL400SR	Reusable Laryngeal Mask Silicone Size 4	1/pack
PMDL500SR	Reusable Laryngeal Mask Silicone Size 5	1/pack



### promed Laryngeal Mask Airways Reusable Reinforced

Item No.	Product	Qty
PMDL100SRR	Reusable Reinforced Laryngeal Mask Silicone Size 1	1/pack
PMDL150SRR	Reusable Reinforced Laryngeal Mask Silicone Size 1.5	1/pack
PMDL200SRR	Reusable Reinforced Laryngeal Mask Silicone Size 2	1/pack
PMDL250SRR	Reusable Reinforced Laryngeal Mask Silicone Size 2.5	1/pack
PMDL300SRR	Reusable Reinforced Laryngeal Mask Silicone Size 3	1/pack
PMDL400SRR	Reusable Reinforced Laryngeal Mask Silicone Size 4	1/pack
PMDL500SRR	Reusable Reinforced Laryngeal Mask Silicone Size 5	1/pack

## promed Intubating Stylets & Oxygen Catheters

The promed intubating stylets are made of a malleable aluminium rod with smooth plastic sheath which can be easily manipulated to the required shape and remains firm when in use. promed oxygen catheters have a round distal eye and eight lateral eyes, a terminal hole and male connector.

Item No.	Product	Qty
PMD07106	Intubating Stylet 06FR pk10	10/pack
PMD07110	Intubating Stylet 10FR pk10	10/pack
PMD07114	Intubating Stylet 14FR pk10	10/pack
PMD06107	Oxygen Catheter Blue 10FG 40cm	100/pack
PMD06108	Oxygen Catheter Blue 14FG 40cm	100/pack



### promed Bacterial/Viral Filter

Item No.	Product	Qty
PMD05000	Bacterial/Viral Filter LuerLock Sampling Prt	50/pack
PMD05001	Bacterial/Viral Filter	50/pack



### promed Combined HME/Bacterial/Viral Filter

PMD05002	HME, Luer-Lock Sampling Port	50/pack
PMD05003	HME with B/V Filter, Luer-Lock Sampling Port	50/pack
PMD05004	HME Filter Luer-Lock Sampling Port Infant	1/pack
PMD05005	HME Filter Luer-Lock Sampling Port Child	1/pack
PMD05006	HME Filter Luer-Lock Sampling Port Adult	1/pack



# promed Respiratory Adaptors & Anaesthesia Accessories

promed Respiratory accessories and adaptors compliment the wide range of respiratory products and allow for customisation of oxygen, gas and medication delivery to meet specific clinical needs.

## promed adaptors

Item No.	Product	Qty
PMD00100	T Adapter	50/pack
PMD00101	Parallel Wye Y	50/pack
PMD00102	Flow Meter Adaptor	50/pack
PMD00103	Tube Connector	250/pack
PMD00104	Universal Wye Connector with Cap	50/pack
PMD00105	Ventilator Y	50/pack
PMD00106	Elbow ACC Standard	50/pack
PMD00107	Acc Aerosol T	50/pack
PMD00108	Adaptor 22mm OD X 15mm ID	50/pack
PMD00109	Oxygen Reducer	50/pack



## promed Anaesthesia Accessories

Item No.	Product	Qty
PMD00001	Elbow with Cap Gas Sampling Luer-Lock	50/Box
PMD00002	Tubing Gas Sampling 3m (10ft)	30/pack
PMD00005	Breathing Bag 1-Liter Neoprene	20/Box
PMD00004	Breathing Bag 2-Liter Neoprene	20/Box
PMD00003	Breathing Bag 3-Liter Neoprene	20/Box



## promed Incentive Spirometer

The promed 3-ball, flow-oriented incentive spirometer features three color-coded balls in chambers to measure inspired volume. The respirometer is used to help improve lung performance in conditions where lung function may be compromised due to various respiratory disorders. This device enables patients to exercise their lungs with the ability to visually examine progress.



- Wide flow rate range from 600-1200mL/sec
- Three colour-coded balls/three chambers
- Minimum flow imprinted on each chamber
- Compact, light-weight design and break-resistant plastic
- Removable tubing for cleaning convenience
- Clear tubing enables clinical inspection of any sectors
- Anatomically designed mouthpiece for patient comfort

Item No.	Product	Qty
PMD00120	Respirometer	20/pack

## promed Yankauer Suction Devices

The promed Yankauer suction handle is comfortable to hold and provides continuous suction. Its strength and durability makes it an ideal choice for airway management suctioning in Emergency, ICU and other ward use. The rigid connector ensures a firm connection to any standard suction.

Item No.	Product	Qty
PMD3404	Yankauer, Vented, Double Wrapped Sterile	50/Box
PMD3405	Yankauer, Non Vented, Double Wrapped Sterile	50/Box

## promed Suction Tubing

The promed Sterile Suction Tubing is designed to be used in conjunction with all suction accessories, including the promed Yankauer Suction Handles. The promed Sterile Suction Tubing is designed to be kink resistant, ensuring suctioning of all fluids including vomitus.

Item No.	Product	Qty
PMD57150	Suction Tubing Sterile, 1.5m Grooved ID 6.5mm OD 10mm	50/Box
PMD57200	Suction Tubing Sterile, 2m Grooved ID 6.5mm OD 10mm	50/Box
PMD57300	Suction Tubing Sterile, 3m Grooved ID 6.5mm OD 10mm	50/Box
PMD57450	Suction Tubing Sterile, 4.5m Grooved ID 6.5mm OD 10mm	50/Box
PMD57600	Suction Tubing Sterile, 6m Grooved ID 6.5mm OD 10mm	50/Box

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

# Order Online NOW!

[www.thermofisher.com.au/godirect](http://www.thermofisher.com.au/godirect)



## goDirect offers:

- 24/7 instant access to over 150,000 products: real-time including YOUR price & stock availability
- No Minimum Order Value (MOV) surcharge for online purchases
- Save time with templates for your shopping lists and frequently purchased products
- Integrated freight tracking with automatic order confirmation, dispatch advice and alerts to estimated time of arrival (ETA)
- Quotes, invoices, statements and reports available online
- Optional authorisation levels for placing and releasing orders for process control



---

**In Australia:**

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, [ordersau@thermofisher.com](mailto:ordersau@thermofisher.com)

**In New Zealand:**

For customer service, call 0800-933-966

To fax an order, use 0800-329-246

To email an order, [ordersnz@thermofisher.com](mailto:ordersnz@thermofisher.com)

Find out more at [thermofisher.com.au/healthcare](https://thermofisher.com.au/healthcare)

**ThermoFisher**  
S C I E N T I F I C