



Radiology Catalogue

- Diagnostic & X-Ray Consumables
- Personal Protective Equipment
- Surgical Consumables
- Ultrasound Consumables

radiology & cath lab
solutions tailored for you

ThermoFisher
SCIENTIFIC

The world leader in serving science

Radiology Solutions for Your Practice

Thermo Fisher Scientific offer a comprehensive range of consumable products designed for a broad range of radiology specialties.

We have long standing partnerships with global market leading suppliers and also complement these key agencies with our own private brand - promed®.

We have refined our vast selection of products into this catalogue to help you provide positive patient outcomes. Some of the speciality products featured in this catalogue include our customised promed® radiology & cath lab procedure packs and drapes as well as aprons and wearables.

X-ray & diagnostic consumables, such as positioning foams and patient transfer mattresses, Parker Laboratories' ultrasound transmission gels, our extensive range of promed® Personal Protective Equipment (PPE) as well as the signosRT® portable ultrasound device for bladder scanning and midwifery applications have also been included. The inclusion of PDI and Briemar infection control consumables complete our portfolio.

Brands we supply



Diagnostic Imaging Solutions

Thermo Fisher Scientific has a long established reputation as a reputable Australian supplier of medical imaging consumable products to the diagnostic imaging and healthcare communities. Our portfolio includes Ultrasound products by Parker Laboratories, Signostics as well as our own brand products for X-Ray consumables such as positioning foams, patient transfer mattresses and personal protective equipment.



Infection Control Solutions

Our flagship brand promed® continues to provide a 'solutions-based' approach to customer needs, specialising in clinical apparel and disinfection. The promed® brand range has a strong reputation for reliable and competitively priced clinical consumables. Complimenting the promed® range, we also offer Briemar® and PDI® branded products.

Our range of infection control products is constantly expanding with new and innovative consumable solutions that compliment the national health policy, professional needs and work practices in the dynamic healthcare community.



Surgical & Operating Room Solutions

The extensive range of promed® surgical drapes, surgical gowns, operating room consumables and procedure packs have been designed in Australia, in consultation with Australian Surgeons, Radiologists, Theatre Nurses and Nurse Unit Managers.

Using a promed® procedure pack for your sterile procedures will:

- Reduce set-up time: fewer products to open, freeing up theatre staff for other duties.
- Streamline ordering & stock control: fewer products to order.
- Reduce Waste: less packaging to dispose.
- Assist maintaining sterility: presented in an aseptic format in line with ACORN guidelines to reduce the risk of contamination.



Ordering online?



Use goDirect...

Go Direct is Thermo Fisher's online portal for the latest product information, stock availability and more, supporting our extensive range of products.



Keyword

Brand

Search

Powerful Search functionality,

Quotes

Easily request quotes

Freight Tracking

Keep track of your order



Some of the best reasons to shop online!

- Powerful Search functionality, to find products easily
- Advanced filtering systems to simplify search results
- 24/7 instant access to over 150,000 products
- Templates for frequently purchased products
- Integrated freight tracking
- Quotes, invoices, statements and reports available online
- Optional authorisation levels for placing and releasing orders

www.thermofisher.com.au/godirect

Contents

Diagnostic Consumables

Biopsy Punches & Grid.....	6
Scalpels & Sharps Safety.....	8
Skin Markers	10

Personal Protective Equipment & Infection Control

Aprons, Light-weight and Lead-free.....	11
Hand Hygiene.....	12
Examination Gloves	13
Sterile Surgical Gowns.....	14
Wipes – patient & equipment.....	15

Patient Positioning & Transfer

Positioning Sponges.....	16
Patient Transfer Mattresses, Sandbags & Sliders	17

Procedure Packs

18

Surgical Consumables

Patient Drapes.....	21
Equipment Drapes.....	24
Surgical Consumables.....	26

Ultrasound Devices

Thermasonic® Gel Warmers.....	30
signosRT® Portable Ultrasound.....	31

Ultrasound Consumables

Ultrasound Gel & Lotion.....	32
Electrode Gel & Skin Preparation	35
Probe Covers – Sterile & Non-Sterile	36

X-Ray Consumables

Image Printing and X-Ray Envelopes	37
X-Ray Markers (Metal Stencils)	38
X-Ray Stickers and Labels.....	40

Diagnostic Consumables

Automatic and Semi-Automatic Biopsy Needles

Gallini is a worldwide leader in developing, manufacturing and marketing healthcare products and has created a wide variety of disposable, high quality biopsy needles for taking a range of tissue specimens within radiology, urology, oncology, gynaecology, incontinence and surgery.

Gallini TZ - Semi-Automatic Biopsy Needle

Features:

- One hand operation
- Light-weight and compact
- Availability of coaxial needle for multiple biopsy


Item Code	Description	UOM
GLNTZ14-10EC	TZ Semi-Automatic Biopsy Needle 14g 10cm	10/Pack
GLNTZ14-16EC	TZ Semi-Automatic Biopsy Needle 14g 16cm	10/Pack
GLNTZ14-20EC	TZ Semi-Automatic Biopsy Needle 14g 20cm	10/Pack
GLNTZ16-10EC	TZ Semi-Automatic Biopsy Needle 16g 10cm	10/Pack
GLNTZ16-16EC	TZ Semi-Automatic Biopsy Needle 16g 16cm	10/Pack
GLNTZ16-20EC	TZ Semi-Automatic Biopsy Needle 16g 20cm	10/Pack
GLNTZ18-10EC	TZ Semi-Automatic Biopsy Needle 18g 10cm	10/Pack
GLNTZ18-16EC	TZ Semi-Automatic Biopsy Needle 18g 16cm	10/Pack
GLNTZ18-20EC	TZ Semi-Automatic Biopsy Needle 18g 20cm	10/Pack
GLNTZ20-10EC	TZ Semi-Automatic Biopsy Needle 20g 10cm	10/Pack
GLNTZ20-16EC	TZ Semi-Automatic Biopsy Needle 20g 16cm	10/Pack
GLNTZ20-20EC	TZ Semi-Automatic Biopsy Needle 20g 20cm	10/Pack

Gallini SPEED-CUT - Automatic Biopsy System

Features:

- Choice of two firing triggers
- Variable throw (15 - 22m)
- One hand operation with two stage arming to reduce effort
- Automatic and semi-automatic procedure
- Availability of coaxial needle for multiple biopsy

Item Code	Description	UOM
GLNSC14-10EC	SPEED-CUT Automatic Biopsy System 14g 10cm	10/Pack
GLNSC14-16EC	SPEED-CUT Automatic Biopsy System 14g 16cm	10/Pack
GLNSC14-20EC	SPEED-CUT Automatic Biopsy System 14g 20cm	10/Pack
GLNSC16-10EC	SPEED-CUT Automatic Biopsy System 16g 10cm	10/Pack
GLNSC16-16EC	SPEED-CUT Automatic Biopsy System 16g 16cm	10/Pack
GLNSC16-20EC	SPEED-CUT Automatic Biopsy System 16g 20cm	10/Pack
GLNSC18-10EC	SPEED-CUT Automatic Biopsy System 18g 10cm	10/Pack
GLNSC18-16EC	SPEED-CUT Automatic Biopsy System 18g 16cm	10/Pack
GLNSC18-20EC	SPEED-CUT Automatic Biopsy System 18g 20cm	10/Pack
GLNSC20-10EC	SPEED-CUT Automatic Biopsy System 20g 10cm	10/Pack
GLNSC20-16EC	SPEED-CUT Automatic Biopsy System 20g 16cm	10/Pack
GLNSC20-20EC	SPEED-CUT Automatic Biopsy System 20g 20cm	10/Pack



Diagnostic Consumables continued

Fast Find Grid

A low-cost positioning grid for quick, accurate CT biopsies and drainage procedures. The Fast Find Grid is a simple disposable tool that provides total accuracy and will find the point of entry FIRST time, EVERY time.

The Fast Find Grid is:

- Accurate
- Quick
- Requires no guesswork
- No repositioning
- Disposable and therefore eliminates cross-contamination

All benefits result in greater comfort for your patient.

SAVE TIME during your interventional procedures by increasing the accuracy of your needle placement cutting down the number of times you need to re-position a needle or re-scan the patient. The Fast Find Grid is useful for all needle biopsies and drainage procedures and can be applied to any part of the patient's anatomy.

Features:

- Grid produces easily identifiable dot pattern on each CT image for positive reference
- Flexible grid conforms to patient's body
- Self-stick tape strips for assured adherence to skin
- Porous material allows marking of patient skin without having to remove the grid 1cm markings for precise needle placement
- Disposable to prevent cross contamination



Item Code	Description	UOM
WBMG6500	Fast Find Grid	25/Pack

Scalpels & Sharp Safety Devices

Canica® Scalpel Handles

The Canica® Standard Scalpel liberates the full level of surgical skill in an ergonomic ejectable blade scalpel handle.

Benefits:

- New level of surgical control and operating room safety
- Push-button blade ejection for minimal disposal risk
- Ergonomically-shaped for enhanced surgical dexterity
- Comes to hand like a fine writing instrument

Item Code	Description	UOM
DRYD6908	Canica® Safety handle size 3 long reusable Non-Sterile	1/Ctn
DRYD6907	Canica® Safety handle size 3 reusable Non-Sterile	1/Ctn
DRYD6910	Canica® Safety handle size 4 long reusable Non-Sterile	1/Ctn
DRYD6905	Canica® Standard scalpel handle rounded reusable	1/Ctn



KAI Scalpel Blades

Benefits:

- Stainless Steel blades
- With transparent cap
- Sterilised by E.O
- 20 scalpels per dispenser box

Item Code	Description	UOM
SMSKAI-NO10	KAI Scalpel blade No.10 Carbon Steel Sterile	100/Pack
SMSKAI-NO11	KAI Scalpel blade No.11 Carbon Steel Sterile	100/Pack
SMSKAI-NO12	KAI Scalpel blade No.12 Carbon Steel Sterile	100/Pack
SMSKAI-NO15	KAI Scalpel blade No.15 Carbon Steel Sterile	100/Pack
SMSKAI-NO20	KAI Scalpel blade No.20 Carbon Steel Sterile	100/Pack
SMSKAI-NO21	KAI Scalpel blade No.21 Carbon Steel Sterile	100/Pack
SMSKAI-NO22	KAI Scalpel blade No.22 Carbon Steel Sterile	100/Pack
SMSKAI-NO23	KAI Scalpel blade No.23 Carbon Steel Sterile	100/Pack
SMSKAI-NO24	KAI Scalpel blade No.24 Carbon Steel Sterile	100/Pack
SMSKAI-NO10SS	KAI Scalpel blade No.10 Stainless Steel, Sterile	100/Pack
SMSKAI-NO11SS	KAI Scalpel blade No.11 Stainless Steel, Sterile	100/Pack
SMSKAI-NO12SS	KAI Scalpel blade No.12 Stainless Steel, Sterile	100/Pack
SMSKAI-NO15SS	KAI Scalpel blade No.15 Stainless Steel, Sterile	100/Pack
SMSKAI-NO20SS	KAI Scalpel blade No.20 Stainless Steel, Sterile	100/Pack
SMSKAI-NO21SS	KAI Scalpel blade No.21 Stainless Steel, Sterile	100/Pack
SMSKAI-NO22SS	KAI Scalpel blade No.22 Stainless Steel, Sterile	100/Pack
SMSKAI-NO23SS	KAI Scalpel blade No.23 Stainless Steel, Sterile	100/Pack
SMSKAI-NO24SS	KAI Scalpel blade No.24 Stainless Steel, Sterile	100/Pack



Scalpels & Sharp Safety Devices continued

Kai Scalpel Handles with Blade

Available in standard or 'safety' model, all feature the sharp quality KAI blades in a range of commonly used sizes.

Item Code	Description	UOM
SMSSCALP10	KAI Scalpel with plastic handle No.10 Sterile	20/Box
SMSSCALP11	KAI Scalpel with plastic handle No.11 Sterile	20/Box
SMSSCALP12	KAI Scalpel with plastic handle No.12 Sterile	20/Box
SMSSCALP15	KAI Scalpel with plastic handle No.15 Sterile	20/Box
SMSSCALP20	KAI Scalpel with plastic handle No.20 Sterile	20/Box
SMSSCALP21	KAI Scalpel with plastic handle No.21 Sterile	20/Box
SMSSCALP22	KAI Scalpel with plastic handle No.22 Sterile	20/Box
SMSSCALP23	KAI Scalpel with plastic handle No.23 Sterile	20/Box
SMSSCALP24	KAI Scalpel with plastic handle No.24 Sterile	20/Box
SMSKAI-810	KAI Safety Scalpel with No. 10 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-811	KAI Safety Scalpel with No. 11 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-812	KAI Safety Scalpel with No. 12 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-815	KAI Safety Scalpel with No. 15 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-820	KAI Safety Scalpel with No. 20 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-822	KAI Safety Scalpel with No. 22 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-823	KAI Safety Scalpel with No. 23 blade Sterile retractable safety scalpel	10/Pack
SMSKAI-824	KAI Safety Scalpel with No. 24 blade Sterile retractable safety scalpel	10/Pack



Sharps Safety

Item Code	Description	UOM
TYC31142576	Devon® 3593 blade shield scalpel holder purple	96/Ctn
PMDBDS001	promed® scalpel transfer tray yellow Sterile	20/Ctn
DRY25-2244	Scrubsafe sharps transfer tray yellow Sterile	24/Ctn
SMSRSUD6495	Sharpak 1.58L	48/Ctn
SMSAMC220-44	Sharpak 220 plus 17L max fill capacity	1/Ctn
SMSRSUD2600	Sharpak 60 plus vol 6L sud 2600	25/Ctn
SMSSUD6201	Sharpak Compak volume 0.5L	24/Ctn
PMDBDS001	promed® scalpel transfer tray yellow Sterile	20/Ctn



SMSRSUD2600



DRY25-2244

Skin Markers

DeRoyal® Pre-Op Mini Markers

These non-sterile mini markers offer convenient storage with every 100 pens coming in their own reusable storage tub. The size and length reduces accidental 'reuse' minimising the opportunity for cross contamination. The ink is non-tattooing and non-toxic.

Item Code	Description	UOM
DRY26-018NS	DeRoyal® Mini-Marker Pre-op skin marker non-sterile 2 tubs of 100pcs non-sterile single use	200/Ctn
DRY26-017NS	DeRoyal® Mini-Marker Pre-op skin marker non-sterile 8 tubs of 100pcs non-sterile single use	800/Ctn



DeRoyal® Sterile Skin Markers

These individually packed sterile skin markers produce smooth, clear line markings on patients skin. Feature prep-resistant gentian violet ink that is non-tattooing and non-toxic.

Item Code	Description	UOM
DRY26-002	DeRoyal® Skin marker and ruler fine tip gentian violet	25/Ctn
DRY26-001	DeRoyal® Skin marker and ruler regular tip gentian violet	25/Ctn



Tyco/Devon™ Sterile Surgical Skin Markers

Item Code	Description	UOM
TYC31145926	Tyco/Devon™ 150 surgical marker regular tip with ruler cap	25/Ctn
TYC31145942	Tyco/Devon™ 151 surgical marker fine tip with ruler cap	25/Ctn
TYC31145900	Tyco/Devon™ 158 surgical marker dual regular and fine tip with ruler cap	25/Ctn
TYC31145819	Tyco/Devon™ 160 surgical marker regular tip gentian violet	25/Ctn
TYC31145785	Tyco/Devon™ 160-R surgical marker regular tip gentian violet with flexi ruler	25/Ctn
TYC31145850	Tyco/Devon™ Duramarker - cs marker red ink	25/Ctn
TYC31145827	Tyco/Devon™ Surgical marker fine tip with flex ruler	25/Ctn

promed® Personal Protective Equipment


Lead Aprons

All staff in a radiographic room during x-ray exposures who are not standing behind protective screens must wear protective clothing. A safety assessment must be undertaken to ensure all staff wear an appropriate lead-equivalent apron giving consideration to the type of procedure and position of the person in relation to the x-ray source.

- Over 12 different designs and styles
- Well respected US manufacturer has been manufacturing lead aprons since 1930
- More than 400 different fabric colours to choose from
- Lead and lead-free options available, depending on budget and personal preference
- Standard and light-weight lead (True Light and Economy range)
- Standard and light-weight lead free (Starlight and Prestige range)

Customised to
suit every shape
and size! 4-6 week
lead time

Item Code	Description	UOM
PMDAU100	promed® Lead protective apron 100cm - unisex	1/Pack
PMDAU090	promed® Lead protective apron 90cm - unisex	1/Pack
PMDAU095	promed® Lead protective apron 95cm - unisex	1/Pack
PMDCU085L	promed® Lead protective comfort apron large straps 85cm - unisex	1/Pack
PMDCU085M	promed® Lead protective comfort apron medium straps 85cm - unisex	1/Pack
PMDCU085S	promed® Lead protective comfort apron small straps 85cm - unisex	1/Pack
PMDCU085XS	promed® Lead protective comfort apron x-small straps 85cm - unisex	1/Pack
PMDCU085XL	promed® Lead protective comfort apron x-large straps 85cm - unisex	1/Pack
PMDUMT0	promed® Lead protective mandarin thyroid cover - unisex	1/Pack
PMDV00LS	promed® Lead protective vest velcro closure small - ladies	1/Pack
PMDV00LM	promed® Lead protective vest velcro closure medium - ladies	1/Pack
PMDV00LL	promed® Lead protective vest velcro closure large - ladies	1/Pack
PMDV00LXL	promed® Lead protective vest velcro closure x-large - ladies	1/Pack
PMDV00SLXXL	promed® Lead protective vest velcro closure xx-large - ladies	1/Pack
PMDV00MS	promed® Lead protective vest velcro closure small - mens	1/Pack
PMDV00MM	promed® Lead protective vest velcro closure medium - mens	1/Pack
PMDV00ML	promed® Lead protective vest velcro closure large - mens	1/Pack
PMDV00MXL	promed® Lead protective vest velcro closure x-large - mens	1/Pack
PMDV00SMXXL	promed® Lead protective vest velcro closure xx-large - mens	1/Pack
PMDS00US	promed® Lead protective skirt velcro closure small - unisex	1/Pack
PMDS00UM	promed® Lead protective skirt velcro closure medium - unisex	1/Pack
PMDS00UL	promed® Lead protective skirt velcro closure large - unisex	1/Pack
PMDS00UXL	promed® Lead protective skirt velcro closure x-large - unisex	1/Pack
PMDS00UXXL	promed® Lead protective skirt velcro closure xx-large - unisex	1/Pack
PMDACUSTOM	promed® Lead protective apron - Customised to order	1/Pack
PMDAPBAG	promed® Lead protective apron carry bag	1/Pack
PMDPPSHIN	promed® Lead protective shin guards with stiffeners	1/Pack
PMDPPEYE	promed® Protective eye wearwear	1/Pack



Still can't find what you are looking for? Call our customer service team on 1300 735 292 to see if we can offer you aprons customised to suit your requirements.

Hand Hygiene

Dermasoft® Sanitising Hand Gels

The range of Dermasoft® Sanitising Hand Gels are soft on skin but hard on bacteria. Dermasoft® Sanitising Hand Gel is not a replacement for hand-washing. It is used in conjunction with a hand washing program to provide the ultimate in hand hygiene. Dermasoft® Sanitising Hand Gel offers a more convenient means of sanitising your hands. As the product requires no water, it can be applied in seconds and eliminates disposal of paper towels.

Independent tests have shown Dermasoft® Sanitising Hand Gel kills 99.99% (Log 4 reduction) of the following organisms in 30 seconds upon application*:

- Staphylococcus aureus
- MRSA (Methicillin Resistant Staphylococcus aureus)
- Pseudomonas aeruginosa
- E. coli
- Salmonella typhimurium
- Candida albicans
- Listeria monocytogenes
- Streptococcus agalactiae

**Microtech Biocidal evaluation of Dermasoft® Sanitising Hand Gel, ref 97036495, 1998*

Item Code	Description	UOM
CCT6782801	Dermasoft® Antibac hand cleanser 1L	8/Pack
CCT16438	Dermasoft® Extra hand wash 1L	8/Pack
CCT16394	Dermasoft® Luxury 1L	8/Pack
CCT17544	Dermasoft® Moisturiser barrier cream 1L	4/Pack
CCT7354991	Dermasoft® Sanitising hand gel 100ml tube	12/Pack
CCT16458	Dermasoft® Sanitising hand gel 1L	4/Pack
CCT17536	Dermasoft® Sanitising hand gel 500ml pump pack	8/Pack



Microshield®

Item Code	Description	UOM
JCMTJJ61222	Microshield® 2 skin cleanser 1.5L	1/Pack
JCMTJJ61354	Microshield® 2 skin cleanser 500ml	1/Pack
JCMTJJ61351	Microshield® 4 hand cleanser/moisturiser 500ml	1/Pack
JCMTJJ61357	Microshield® Handrub 500mL	1/Pack
JCMTJJ61224	Microshield® Handwash 1.5L	1/Pack

promed® Examination Gloves

Item Code	Description	UOM
PMD173060	promed® Nitrile Soft Examination Gloves; X-small	200/Box
PMD173070	promed® Nitrile Soft Examination Gloves; Small	200/Box
PMD173080	promed® Nitrile Soft Examination Gloves; Medium	200/Box
PMD173090	promed® Nitrile Soft Examination Gloves; Large	200/Box
PMD173095	promed® Nitrile Soft Examination Gloves; X-large	200/Box
PMD172060	promed® Nitrile Long Cuff Examination Gloves; X-small	100/Box
PMD172070	promed® Nitrile Long Cuff Examination Gloves; Small	100/Box
PMD172080	promed® Nitrile Long Cuff Examination Gloves; Medium	100/Box
PMD172090	promed® Nitrile Long Cuff Examination Gloves; Large	100/Box
PMD172095	promed® Nitrile Long Cuff Examination Gloves; X-large	100/Box
PMD148060	promed® Latex Powder Free Examination Gloves; X-small	100/Box
PMD148070	promed® Latex Powder Free Examination Gloves; Small	100/Box
PMD148080	promed® Latex Powder Free Examination Gloves; Medium	100/Box
PMD148090	promed® Latex Powder Free Examination Gloves; Large	100/Box
PMD148095	promed® Latex Powder Free Examination Gloves; X-large	100/Box
PMD1300	promed® Vinyl Powder Free Examination Gloves; Small	100/Box
PMD1301	promed® Vinyl Powder Free Examination Gloves; Medium	100/Box
PMD1302	promed® Vinyl Powder Free Examination Gloves; Large	100/Box
PMD1303	promed® Vinyl Powder Free Examination Gloves; X-large	100/Box
APD3210	promed® Glove Holder, Wall Mounted One Tier	EA
APD3450	promed® Glove Holder, Wall Mounted Two Tier	EA
APD3640	promed® Glove Holder, Wall Mounted Three Tier	EA
APD3740	promed® Glove Holder, Wall Mounted; For Nitrile Long Cuff Three Tier	EA



promed® Surgical Gowns

promed® Standard Surgical Gowns

The promed® Standard Surgical Gown is made with 50gsm SMMS to provide maximum fluid repellence and increased wearer comfort. With ultrasonically welded sleeve seams it meets AAMI Level 3 standards making it ideal for use across a variety of areas and procedures.

All promed® Standard Surgical Gowns come with a generous outer wrap, and absorbent cotton huck towel.



Item Code	Description	Size	UOM
PMD8221102	promed® Standard Surgical Gowns	Small	25/Ctn
PMD8231102	promed® Standard Surgical Gowns	Medium	25/Ctn
PMD8241102	promed® Standard Surgical Gowns	Large	25/Ctn
PMD8251102	promed® Standard Surgical Gowns	X-Large	25/Ctn
PMD8261102	promed® Standard Surgical Gowns	XX-Large	25/Ctn
PMD8271102	promed® Standard Surgical Gowns	Large X-Long	25/Ctn
PMD8281102	promed® Standard Surgical Gowns	X-Large X-Long	25/Ctn

promed® Impervious Reinforced Surgical Gowns

The promed® Impervious Reinforced Surgical Gown is made with a highly breathable, and comfortable 68gsm blue spunlace. The reinforcement material is a 20gsm polypropylene breathable film which combined with the ultrasonically welded sleeves improves wearer comfort and safety by reducing the potential for fluid strikethrough.

All promed® Impervious Reinforced Surgical Gowns come with a generous outer wrap, and absorbent cotton huck towel.



Item Code	Description	Size	UOM
PMD8421102	promed® Impervious Reinforced	Small	25/Ctn
PMD8431102	promed® Impervious Reinforced	Medium	25/Ctn
PMD8441102	promed® Impervious Reinforced	Large	25/Ctn
PMD8451102	promed® Impervious Reinforced	X-Large	25/Ctn
PMD8461102	promed® Impervious Reinforced	XX-Large	25/Ctn
PMD8471102	promed® Impervious Reinforced	Large X-Long	25/Ctn
PMD8481102	promed® Impervious Reinforced	X-Large X-Long	25/Ctn

Wipes

PDI® Hygea Patient Wipes

A hand and face wipe specially designed to keep patients clean and refreshed.

- Hypoallergenic and dermatologically tested, paraben and alcohol free
- Label offers space for patient/ward details
- Significantly reduces bathing time, laundry cost and staff time
- Deep-well lid to avoid cross-contamination
- Wall brackets available



Item Code	Description	UOM
PDIXP00195	PDI® Hygea Bodycare 60 count flat pack	8/Ctn
PDIXP00162	PDI® Hygea Hands & face 125 count canister	6/Ctn
PDIXP00161C	PDI® Hygea Hands & face 200 count canister	6/Ctn
PDIXP00164	PDI® Hygea Hands & face 24 count flat pack	24/Ctn
PDIXP00166	PDI® Hygea Hands & face 80 count dispenser box	6/Ctn
PDIXP00165	PDI® Hygea Hands & face 80 count refill	12/Ctn
PDIXP00207	PDI® Hygea Moist toilet tissue 10 count pocket	24/Ctn
PDIXP00208	PDI® Hygea Refreshing wipe 10 count pocket	24/Ctn
PDIXP00287	PDI® Hygea Newborn babywipes 125 count canister	6/Ctn
PDIXP216	PDI® Hygea Bodycare 46g 100 count dry wipes	12/Ctn
PDIXP217	PDI® Hygea Bodycare 60g 70 count dry wipes	12/Ctn
PDIXP00299-P	PDI® Hygea Purple bracket for Hygea Patient Care Canister Wipes	1/Box

Hard Surface Wipes

Our wipe range is an efficient and cost effective method for the cleaning of medical devices and surfaces within the healthcare environment. Each wipe offers low lint, no messy mixing or rinsing and is competitively priced.

promed® Multi-Purpose wipes

Item Code	Description	UOM
PMDNW7100	promed® Protowel disposable towel 33×32cm white	1200/Pack
PMDNW7000	promed® Prowipe - antibacterial wipes 75 count	12/Pack
PMDNW8000	promed® Neutrawipe - neutral detergent wipe 75s	12/Ctn
PMDNW8001	promed® Neutrawipe - neutral detergent wipe 140s	4/Ctn



PDI® Sani-Cloth Detergent Wipes

















Item Code	Description	UOM
PDIXP00230	PDI® Sani-cloth detergent 100 count flat pack	6/Ctn
PDIXP00301C	PDI® Sani-cloth detergent 150 count canister	6/Ctn
PDIXP00175	PDI® Sani-cloth detergent 225 count bucket	1/Box
PDIXP00210	PDI® Sani-cloth detergent 225 count refill	4/Ctn
PDIXP00186	PDI® Sani-cloth detergent 60 count flat pack	8/Ctn
PDIXP00299-Y	PDI® Yellow bracket for Sani-Cloth Detergent Wipes	1/Box




Patient Positioning

Positioning Sponges

Proper patient positioning is essential in order to obtain accurate diagnostic images for your patient. Thermo Fisher offer several solutions to facilitate you obtaining films of premium quality that aid diagnosis. Our positioning sponges are made of high density foam in both uncovered or easy to clean PVC covered options. Various sizes are available including custom cut options:

Item Code	Description		UOM
JCBFP13-FV001390	Stepped Finger 130mm × 150mm × 60mm		EA
JCBFP10C (PVC covered) JCBFP10-FV001265	Large Triangle 457mm × 159mm × 159mm		EA
JCBFP11C (PVC covered) JCBFP11-FV001565	Triangle 254mm × 127mm × 127mm		EA
JCBFP14C (PVC covered) JCBFP14-FV001389	Small Triangle 140mm × 170mm × 100mm		EA
JCBFP15C (PVC covered) JCBFP15-FV001514	Channel 280mm × 152mm × 152mm		EA
JCBFP17C (PVC covered) JCBFP17-FV001388	Jumbo CT 300mm × 450mm × 175mm		EA
JCBFP1C (PVC covered) JCBFP1-FV001264	Circle 152mm diameter × 51mm		EA
JCBFP2C (PVC covered) JCBFP2-FV001669	Small Triangle 178mm × 102mm × 102mm		EA
JCBFP3C (PVC covered) JCBFP3-FV001267	Large 45° Triangle 254mm × 178mm × 178mm		EA
JCBFP3LC (PVC covered) JCBFP3L-FV001268	45° Triangle Long 460mm × 178mm × 178mm		EA
JCBFP4C (PVC covered) JCBFP4-FV001563	Large Triangle 254mm × 178mm × 178mm		EA
JCBFP5C (PVC covered) JCBFP5-FV001269	Wedge 20° 254mm × 178mm × 75mm		EA
JCBFP6C (PVC covered) JCBFP6-FV001270	Quad 254mm × 178mm × 50mm		EA
JCBFP7C (PVC covered) JCBFP7-FV001670	Cylinder 178mm diameter × 305mm		EA
JCBFP8C (PVC covered) JCBFP8-FV019410	Skull 229mm × 178mm × 152mm		EA
JCBFP9C (PVC covered) JCBFP9-FV001516	Large Block 264mm × 127mm × 127mm		EA



Patient Positioning continued

Further Patient Positioning options available include

Item Code	Description	UOM
JCBFPTU-FV009906	T-vaginal 400mm × 450mm × 100mm (for endocavity ultrasound)	1/Pack
JCBFPFV906073	15 Degree Wedge foam pad 150 × 370 × 100(380)mm	1/Pack
JCBFPFV906073C	15 Degree Wedge foam pad 150 × 370 × 100(380)mm Covered	1/Pack
JCBFPFV906074	25 Degree Wedge foam pad 150 × 370 × 170(410)mm	1/Pack
JCBFPFV906074C	25 Degree Wedge foam pad 150 × 370 × 170(410)mm Covered	1/Pack
JCBFPFV906733	Foam chest pad 430 × 350 × 200mm	1/Pack
JCBFPFV906733C	Foam chest pad 430 × 350 × 200mm Covered	1/Pack
JCBFPFV906772	Foam arm board 600 × 155 × 200mm	1/Pack

Still can't find what you are looking for? Call our customer service team on 1300 735 292 to see if we can have a positioning sponge custom cut to your requirements.

(Please note that the custom cut products incur additional charges)

Patient Transfer Mattresses

Item Code	Description	UOM
JCBMS10-CUSTOM	Custom foam mattress 2100 × 600 × 25mm	1/Pack
JCBFPFV803447C	Foam mattress 1800 × 480 × 20mm covered	1/Pack
JCBFPFV180060040C	Foam mattress 1800 × 600 × 40mm covered heavy duty PVC	1/Pack
JCBFPFV905234C	Foam mattress 1800 × 600 × 50mm covered heavy duty PVC	1/Pack
JCBFPFV805123C	Foam mattress 2100 × 800 × 50mm covered	1/Pack
JCBFPFV220605C	Foam mattress 2200 × 600 × 50mm covered H/D	1/Pack
JCBFPFV220705C	Foam mattress 2200 × 700 × 50mm covered H/D	1/Pack
JCBFPFV044055C	Foam mattress 2200 × 640 × 50mm covered heavy duty PVC	1/Pack
JCBMPS10GNC	X-ray mattress 1800 × 600 × 22mm covered with green nylon fabric	1/Pack
JCBMPS10BGNC	X-ray mattress 1800 × 600 × 22mm with belts covered with green nylon fabric	1/Pack
JCBMPS10	X-ray mattress 1800 × 600 × 25mm	1/Pack
JCBMPS50	X-ray mattress 1800 × 600 × 50mm	1/Pack
JCBMTTM	Trauma transfer mattress 1800 × 500 × 50mm	1/Pack
JCBMTTML	Trauma transfer mattress 1800 × 750 × 50mm	1/Pack

Sandbags

Item Code	Description	UOM
JCBMSB10	1kg sandbag 14cm × 24cm finished size	1/Pack
JCBMSB30	3kg sandbag 19cm × 34cm finished size	1/Pack
JCBMSB60	6kg sandbag 30 × 35cm finished size	1/Pack

Sliders

Item Code	Description	UOM
JCBMSLIDER24X30P	Slider 24 × 30cm standard configuration	1/Pack
JCBMSLIDER35X43L	Slider 35 × 43cm horizontal or landscape black	1/Pack
JCBMSLIDER35X43V	Slider 35 × 43cm vertical or portrait configuration	1/Pack

promed® Procedure Packs

promed® procedure packs

Using a promed® procedure pack in your Radiology practice or Cath lab is the solution to:

- **Reduce set-up time:** fewer products to open, freeing up theatre staff for other duties
- **Streamline ordering & stock control:** fewer products to order
- **Reduce waste:** less packaging to dispose of
- **Cost control efficiency:** minimise the number of component costs per procedure to track
- **Maintaining sterility:** presented in an aseptic format in line with ACORN guidelines to reduce the risk of contamination
- **Traceability:** all components included in a promed® procedure pack can be track back to a specific manufacturer batch



promed® procedure packs can be customised to suit based on surgeon preference, and individual healthcare facility protocol. The following packs have been pre-designed with other leading healthcare providers and may offer a total solution, or be the ideal basis for commencing the design of your own customised promed® procedure pack solution.

promed® Angiography Kit

Item Code: PMD6912

UOM: 4/Ctn

Contains:

4 × Absorbent Wipe / Hand Towel (32.5cm × 60cm)	10 × Gauze Swabs 10cm × 10cm 8ply (Pkt5)
1 × Surgical Scalpel Blade & Handle # 11	2 × Needle 19g × 38mm
1 × Needle 23g × 32mm	1 × Needle Arterial 18g
3 × Prep Sponge	2 × Angiography Drug Label Sheet
1 × IV Set VQS MLL 150cm	1 × Guidewire 145cm 3mm
1 × IV Set FC NVS	1 × Compensator Manifold
1 × Pressure Monitoring Line 60cm	1 × HP Line 120cm 1000PSI
1 × Control Syringe	4 × Syringe Luer Slip 10ml
1 × Gallipot 150ml	1 × Kidney Dish 700ml Clear
1 × Bowl 500ml	1 × Tray Blue
1 × Wrap 100cm × 150cm SMS	1 × I. I. Cover
1 × Lead Screen Cover	1 × Angiography Drape
1 × Backtable Cover / Outer Wrap	

promed® Angiography Add On Kit

Item Code: PMD6948

UOM: 15/Ctn

Contains:

1 × Clamp Forceps 11cm	1 × Drawing Up Needle Blunt 18g × 38mm
1 × Pressure Monitoring Line 60cm	2 × Stopcock 3 Way Leur Lock
1 × Syringe Leur Lock 20ml	1 × Syringe Leur Slip 2ml
1 × Syringe Leur Slip 10ml	1 × Inflation Device
1 × Wrap SMS 100cm × 150cm	1 × Backtable Cover / Outer Wrap

promed® Procedure Packs continued

promed® Electrophysiology Pack

Item Code: PMD6949

UOM: 6/Ctn

Contains:

4 × Blue Cotton Huck / Hand Towel	2 × Drawing Up Needle Blunt 18g × 38mm
2 × Needle 23g × 32mm	1 × Needle Artery 18g × 7cm
4 × Syringe Leur Slip 10ml	1 × Surgical Scalpel Blade & Handle # 11
4 × Clamp Forceps 11cm	3 × Prep Swab / Sponge
10 × Gauze Swabs 10cm × 10cm 8ply (Pkt5)	1 × Drug Label 'Heparin'
2 × Drug Label 'Lignocaine'	1 × Wrap SMS 100cm × 150cm
1 × I.I Cover 60cm × 60cm	1 × I.I Cover 98cm × 86cm
1 × Kidney Dish Clear 700ml	1 × Gallipot Clear 150ml
1 × 4-Compartment Tray	1 × Brachial Drape 75cm × 100cm
1 × Angio Drape with Square Fenestrations	1 × Backtable Cover / Outer Wrap



promed® Electrophysiology Add On Kit

Item Code: PMD6950

UOM: 10/Ctn

Contains:

1 × Giving Set Leur Lock	1 × Pressure Monitoring Line M/F 120cm
1 × Giving Set	1 × Control Syringe 10ml
1 × Compensator Manifold	1 × Wrape SMS 90cm × 90cm



promed® Pacemaker Instrument Pack

Item Code: PMD6952

UOM: 10/Ctn

Contains:

2 × Prep Swab Stick	1 × Forcep Adson Toothed
2 × Forcep Artery Mosquito Curved 12.5cm	2 × Criles Straight 14cm
1 × Forcep Spencer Wells Curved	1 × Gauze Swabs 10cm × 10cm 8ply (Pkt5)
1 × Ligature Needle	1 × Metzenbaum Scissor Curved 16cm
1 × Needle Holder	2 × Forceps Plastic
1 × Scissor Iris Sharp/Sharp 11.5cm	1 × Towel Clamp
1 × Mollison Retractors Sharp 16cm	1 × Kidney Dish Clear 700ml
2 × Gallipot Clear 150ml	1 × Forcep Dressing
1 × 3-Compartment Tray	1 × Bowl 1L
1 × Wrap Heavy Duty SMS 130cm × 150cm	

promed® Procedure Packs continued

promed® PICC Pack Radiology

This multi bundle kit is ideal for handing off relevant component to multiple specialties. Also includes a sterile probe cover.

Item Code: PMD6954

UOM: 5/Ctn

Bundle 1 - Radiology

Contains:

4 × Blue Cotton Huck / Hand Towel	20 × Gauze Plain 10cm × 10cm (Pkt5)
1 × Forcep Spencer Wells Straight 14cm	1 × Drug Label 'Contrast' Yellow
1 × Kidney Dish Clear 700ml	4 × Gallipot Clear 150ml
1 × Wrap SMS 120cm × 120cm	

Bundle 2 - Angio

Contains:

1 × Blue Cotton Huck / Hand Towel	1 × Scissors Sharp/Sharp 13cm
1 × Scissors Sharp/Blunt 12.5cm	1 × Forcep Artery Mosquito Curved 12.5cm
2 × Dressing Tray	1 × Wrap SMS 90cm × 90cm

Bundle 3

Contains:

1 × Brachial Drape 75cm × 100cm	1 × Aquasonic Gel Sachet
1 × Tegaderm Dressing 10cm × 12cm	1 × Syringe Leur Slip 10ml
2 × Syringe Leur Lock 10ml	1 × Needle 23g × 32mm
3 × Drawing Up Needle Blunt 18g × 38mm	2 × Micro Clave Connector
1 × Rubber Band	1 × Wrap SMS 120cm × 120cm



promed® Radiology Pack

Item Code: PMD6955

UOM: 15/Ctn

Contains:

1 × Blue Cotton Huck / Hand Towel	1 × Gauze Plain 10cm × 10cm (Pkt5)
1 × Forcep Spencer Wells Straight 13cm	1 × Drug Label 'Contrast'
4 × Gallipot Clear 150ml	1 × Kidney Dish Clear 700ml
1 × Wrap SMS 120cm × 120cm	



promed® Lumbar Puncture Pack

Item Code: PMD6956

UOM: 15/Ctn

Contains:

1 × Terri Wipe	1 × Gauze Sponge 5cm × 5cm (Pkt5)
1 × Forceps Green	1 × Needle 25g × 16mm
1 × Syringe Leur Lock 5ml	1 × Spinal Manometer
1 × IV3000 Dressing 6cm × 8.5cm	1 × Drug Label 'Contrast'
1 × Specimen Jar with Label 5ml	1 × Wrap SMS 90cm × 90cm
4 × Specimen Jar 10ml with V Bottom (numbered)	



promed® Drapes

promed® Patient Drapes

The extensive range of promed® Surgical Drapes have been designed in Australia, in consultation with Australian Surgeons, Nurses and Theatre Managers.

From the generous sizing of the Promed Drapes, to the absorbency and strength of the materials used, the entire range is designed to meet Australian Theatre Protocols, Standards and Operating Procedures.


The latest quality materials are combined in the manufacture of the Promed Drapes to ensure:

- Absorbency and Fluid Management
- Waterproofing where required to eliminate strikethrough
- Patient comfort through low weight fabrics and skin friendly adhesives
- Operating Theatre safety through low flammability
- Strength and drapability



promed® Angiography Drapes

Item Code	Description	UOM
PMD7150102	promed® Angio/pacemaker drape 370cm × 210cm dual fenestration Non-Sterile	10/Ctn
PMD5329.2	promed® Angiography Drape 180cm × 250cm Sterile	20/Ctn
PMD79711	promed® Angiography Drape with windows 221cm × 343cm Non-Sterile	18/Ctn
PMD7201102	promed® Angiography Drape with comfasorb 210cm × 317cm Sterile	10/Ctn
PMD7200102	promed® Angiography Drape ProDrape 195cm × 300cm Sterile	10/Ctn




PMD7200102

promed® Pacemaker Drapes

A choice of individual drapes, or save on your set up time with our pacemaker drape pack.

Item Code	Description	UOM
PMD7100102	promed® Pacemaker drape ProDrape dual fenestrations 180cm × 320cm Sterile	10/Ctn
PMD5294	promed® Permanent Pacing Drape Sterile	10/Ctn



PMD7201102

promed® Pacemaker Drape Pack

Item Code: PMD6951	UOM: 10/Ctn
Contains:	
1 × Cotton Hand Towel Blue	1 × Angio Drape
1 × I.I Cover 60cm × 60cm	1 × Backtable Cover / Outer Wrap
1 × Wrap SMS 100cm × 150cm	

Minor Procedure Drapes

All feature adhesive around fenestration.

Medidrape


Medidrape is an absorbent two layer material that is impervious to fluids.

Item Code	Material	Drape Size	Fenestration Size	UOM
PMD5302	Medidrape	50cm × 60cm	6.5cm circular	50/Ctn
PMD5662	Medidrape	50cm × 60cm	10cm circular	50/Ctn
PMD5704	Medidrape	60cm × 60cm	4cm × 3cm	50/Ctn
PMD5706	Medidrape	90cm × 90cm	8cm × 4cm	50/Ctn
PMD5705	Medidrape	90cm × 90cm	10cm × 4cm	50/Ctn
PMD5707	Medidrape	90cm × 90cm	10cm × 25cm	50/Ctn
PMD5126S	Medidrape	110cm × 150cm	10cm × 15cm	50/Ctn

ProDrape

ProDrape is impervious to fluids, with a bottom layer designed for patient comfort against the skin.

Item Code	Material	Drape Size	Fenestration Size	UOM
PMD1220102	ProDrape	50cm × 60cm	6.5cm round	50/Ctn
PMD1221102	ProDrape	80cm × 80cm	7cm × 17cm	25/Ctn
PMD1401102	ProDrape	60cm × 90cm	6cm circular	50/Ctn
PMD1234102	ProDrape	60cm × 90cm	6cm circular	50/Ctn
PMD1402102	ProDrape	90cm × 90cm	6cm circular	25/Ctn
PMD1403102	ProDrape	90cm × 120cm	6cm perforated circular	50/Ctn
PMD1236102	ProDrape	90cm × 120cm	6cm perforated circular	50/Ctn
PMD1235102	ProDrape	90cm × 90cm	6cm perforated circular	50/Ctn
PMD1231102	ProDrape	120cm × 120cm	10cm × 15cm	25/Ctn
PMD1404102	ProDrape	120cm × 120cm	6cm circular	50/Ctn
PMD1225102	ProDrape	140cm × 180cm	7cm × 16cm	25/Ctn
PMD1226102	ProDrape	150cm × 80cm	5cm × 15cm	25/Ctn
PMD1222102	ProDrape	152cm × 110cm	10cm × 15cm	25/Ctn
PMD1500102	ProDrape	170cm × 270cm	oval adhesive	10/Ctn




PMD1402102

SMS

SMS provides strength & fluid repellency, with the additional benefit of breathability for patient comfort.

Item Code	Material	Drape Size	Fenestration Size	UOM
PMD5293.0	SMS	40cm × 40cm	7cm circular	50/Ctn
PMD1320102	SMS	60cm × 60cm	6.5cm circular	50/Ctn
PMD5610	SMS	60cm × 90cm	5cm × 15cm	50/Ctn
PMD5611	SMS	80cm × 80cm	10cm × 10cm	50/Ctn
PMD1310102	SMS Heavy Duty	70cm × 90cm	15cm × 5cm	50/Ctn
PMD1312102	SMS Heavy Duty	110cm × 130cm	5cm × 15cm	50/Ctn
PMD1311102	SMS Heavy Duty	120cm × 75cm	15cm × 5cm	50/Ctn



PMD5610

Minor Procedure Drapes continued

promed® Clear Fenestrated Drapes

Made from a strong and clear impervious material these promed® Clear Fenestrated Drapes are ideal for use in a number of procedures where a large visual radius of the surgical site is required.

All drapes are individually packaged in a sterile easy-to-open peel pouch.

Item Code	Drape Size	Fenestration Size	UOM
PMD1126102	50cm × 124cm	10cm × 15cm	50/Ctn
PMD1120102	84cm × 90cm	6cm × 10cm	50/Ctn
PMD1121102	85cm × 136cm	8cm × 21cm	50/Ctn
PMD1124102	100cm × 100cm	5cm × 5cm	50/Ctn
PMD1125102	120cm × 120cm	10cm × 15cm	50/Ctn
PMD1123102	124cm × 180cm	10cm × 15cm	25/Ctn

promed® Bed & Draw Sheets

These promed® Medidrape is ideal for under patient use. They are absorbent, anti-slip and impervious to fluid.

Item Code	Description	UOM
PMD2073	promed® Drape Medidrape 60cm × 100cm blue Non-Sterile	100/Ctn
PMD2072	promed® Drape Medidrape 60cm × 100cm green Non-Sterile	100/Ctn
PMD5049	promed® Drape Medidrape 150cm × 100cm Sterile	50/Ctn
PMD5406	promed® Drape Medidrape 150cm × 180cm Sterile	50/Ctn
PMD5408	promed® Drape Medidrape 180cm × 250cm Sterile	10/Ctn
JCBHS02002122	promed® Waterproof Bed Sheet Blue 1200mm × 50mt	1/Pack

promed® Equipment Drapes

promed® Equipment Drapes

Item Code	Description	UOM
PMD5726BB	promed® B Braun II Cover Non-Sterile 65cm × 80cm	100/Ctn
PMD5731BB	promed® B Braun II Cover Non-Sterile 91cm × 91cm Non-Sterile	60/Ctn
PMD5733BB	promed® B Braun Light Handle Cover Non-Sterile 15cm × 26cm	50/Ctn
PMD29900	promed® Cable Cord Cover	100/Ctn
PMD8220102	promed® Camera sleeve closed end clear 15cm × 200cm	20/Ctn
PMD5483	promed® Drill sleeve 8cm × 100cm	50/Ctn
PMD5347.0	promed® Drill sleeve dexter green	50/Ctn
PMD5122NS	promed® Front Panel Drape Clear Plastic with Double Sided Adhesive Strip Non-Sterile 56cm × 46cm	50/Ctn
PMD5832NS	promed® Burns Chair Cover and Absorbent Panel Non-Sterile 80cm × 180cm	10/Ctn
PMD5485	promed® Drill Sleeve Sterile 8cm × 100cm	50/Ctn
PMD5423	promed® Drill Sleeve Sterile 9cm × 92cm	50/Ctn
PMD5407	promed® Ice Machine Cover with circular fenestration Sterile 150cm × 180cm	25/Ctn
PMD1229	promed® Ice pack cover with flap Sterile	500/Ctn
PMD8222102	promed® II camera / image intensifier sleeve clear poly Sterile 15cm × 132cm	50/Ctn
PMD8202102	promed® II arm & head cover w/tabs 75cm (w) × 140cm (h)	30/Ctn
PMD8203102	promed® Image cover cap bouffant 50cm	40/Ctn
PMD8205102	promed® Image cover cap 75cm × 100cm bouffant	40/Ctn
PMD8206102	promed® Image deep cover cap bouffant 115cm × 135cm	40/Ctn
PMD8204102	promed® Image deep cover cap bouffant 65cm × 80cm	40/Ctn
PMD8201102	promed® Image intensifier arm & head cover Sterile 75cm (w) × 140cm (h)	30/Ctn
PMD5414.0	promed® Image intensifier cover bouffant clear plastic Sterile 50cm	100/Ctn
PMD8215102	promed® Image intensifier ceiling cover Sterile 105cm (w) × 105cm (h)	40/Ctn
PMD8207102	promed® Image intensifier ceiling cover Sterile 109 (w) × 190cm (h)	40/Ctn
PMD5209	promed® Image Intensifier Cover 75cm diameter Sterile	100/Ctn
PMD5216	promed® Image Intensifier Cover 80cm diameter Sterile	100/Ctn
PMD5583NS	promed® Image Intensifier Cover Non-Sterile 90cm × 106cm	50/Ctn
PMDIISC	promed® Image Intensifier Cover - Shower Cap Style	50/Ctn
PMD8208102	promed® Image intensifier screen cover Sterile 60cm (w) × 60cm (h)	40/Ctn
PMD8209102	promed® Image intensifier screen cover Sterile 75cm (w) × 50cm (h)	40/Ctn
PMD8210102	promed® Image intensifier screen cover Sterile 75cm (w) × 80cm (h)	40/Ctn
PMD8213102	promed® Image intensifier screen cover Sterile 95cm (w) × 50cm (h)	40/Ctn
PMD8214102	promed® Image intensifier screen cover Sterile 95cm (w) × 75cm (h)	40/Ctn
PMD8216102	promed® IV stand cover Sterile 23cm × 110cm	40/Ctn
PMD5467	promed® Laser Sleeve Sterile 25cm × 75cm	50/Ctn
PMD5416	promed® Lead screen 70cm × 100cm Sterile	50/Ctn
PMD8218102	promed® Lead screen cover 104cm × 150cm Sterile	20/Ctn
PMD932	promed® Microscope Cover for Zeiss MD Sterile 103cm × 259cm	10/Ctn
PMDMK974	promed® Microscope Cover for Zeiss 65mm (22in) Dome	30/Ctn
PMD5488	promed® Microscope Covers (pair) open ended cuff 20cm × 15cm	50/Ctn



II Cover range



promed® Equipment Drapes continued

promed® Equipment Drapes continued

Item Code	Description	UOM
PMD9009	promed® Orthopaedic Navigation CR2 Batteries Pack - each pack contains 3 batteries	20/Ctn
PMD5241.0	promed® Plastic equipment cover Sterile 104cm × 165cm	25/Ctn
PMD5702	promed® Universal C-Arm Cover Sterile	25/Ctn
PMD5459	promed® Universal Sleeve Covers Sterile (Pairs)	30/Ctn
PMDC337	promed® X-Ray Cassette Bag Sterile	50/Ctn
PMDC637NS	promed® X-Ray Cassette Bag Non-Sterile	50/Ctn



promed® Backtable / Trolley Covers

A comprehensive range provides an ideal solution to suit any sized trolley. Custom sizes are also available upon request.

Item Code	Description	UOM
PMD8111102	promed® Back Table / Trolley Cover Clear poly 122cm × 137cm Sterile	25/Ctn
PMDMC5419F	promed® Back Table / Trolley Cover 137cm × 190cm Sterile	25/Ctn
PMD5419IC	promed® Back Table / Trolley Cover 152.5cm × 228.5cm	25/Ctn
PMD5419S	promed® Back Table / Trolley Cover blue poly 140cm × 190cm Sterile	25/Ctn
PMD5417	promed® Back Table / Trolley Cover blue poly reinforced 190cm × 109cm Sterile	25/Ctn
PMD5419I	promed® Back Table / Trolley Cover heavy duty 137cm × 190cm Sterile	25/Ctn
PMD8103102	promed® Back Table / Trolley Cover heavy duty 137cm × 190cm with absorbent panel Sterile	25/Ctn
PMD8110102	promed® Back Table / Trolley Cover heavy duty 137cm × 227cm with absorbent panel Sterile	20/Ctn
PMD5419HDNS	promed® Back Table / Trolley Cover heavy duty 150cm × 190cm Non-Sterile	75/Ctn
PMD5372.0	promed® Back Table / Trolley Cover heavy duty comfasorb 185cm × 55cm Sterile	25/Ctn



PMD5419S



PMD5417

promed® Clear Utility Drapes

These clear sterile drapes are impervious and available in a range of sizes. All drapes are individually packaged, in a sterile easy-to-open peel pouch.

Item Code	Drape Size	UOM
PMD1102102	60cm × 45cm	50/Ctn
PMD1103102	61cm × 91cm	50/Ctn
PMD1107102	61cm × 183cm	50/Ctn
PMD1104102	76cm × 152cm	50/Ctn
PMD3008	80cm × 90cm	50/Ctn
PMD1105102	90cm × 120cm	50/Ctn
PMD1101102	100cm × 100cm	50/Ctn
PMD3006	150cm × 150cm	50/Ctn

Surgical Consumables

Needle Counters

Our extensive range of Needle Counters contains a number of safety features to assist in preventing a sharps injury. Available in a range of styles and colours including magnetic, foam or a combination.

Item Code	Description	UOM
PMD5026M	promed® Mini Needle Counter Magnet and Foam Sterile	50/Ctn
PMD5010	promed® Needle Counter 10 count Sterile	50/Ctn
PMD5025	promed® Needle Counter 25 Count Sterile	50/Ctn
PMD5005	promed® Needle Counter double Magnet No No's Sterile	50/Ctn
PMD5002	promed® Needle Counter Foam Sterile	72/Ctn
PMD61225S	promed® Needle Counter Magnet & Foam	50/Ctn
PMD4678I	promed® Needle Counter Sterile	40/Ctn
PMD5005M	promed® Needle Counter Mini Double Magnet	50/Ctn
TYC31142410	Devon® 1100 needlecounter single magnet red with clear top 10 count	21/Ctn
TYC31142212	Devon® 1105 Needlecounter 15 ct foam strip Clear top	12/Ctn
TYC31142196	Devon® 1200 needlecounter double magnet red 20 count	21/Ctn
TYC50091201	Devon® 1201 needlecounter double magnet yellow without numbers	21/Ctn
TYC31142147	Devon® 1230 needlecounter double magnet red 30 count	12/Ctn
TYC50091230	Devon® 1230 needlecounter double magnet yellow 30 count	12/Ctn
TYC50091260	Devon® 1260 needlecounter double magnet yellow 60 count	8/Ctn
TYC31142303	Devon® 1330 needlecounter double foam strip red 30 count	96/Ctn
TYC31142238	Devon® 1540 needlecounter bladeguard double white mag red 30 count	12/Ctn
TYC50091614	Devon® 1614 needlecounter bladeguard foam magnet strips yellow 15 count	12/Ctn
TYC31142295	Devon® 1614 needlecounter bladeguard foam/magnet strip red 30 count	12/Ctn
TYC31142451	Devon® 1630 needlecounter bladeguard double foam strips red 30 count	12/Ctn
TYC50091822	Devon® 1822 needlecounter foam adhesive yellow 20 count	12/Ctn
TYC31142402	Devon® 1842 needlecounter foam adhesive 40 count/70 capacity red	8/Ctn
TYC31142394	Devon® 1910 needlecounter bladeguard single foam magnet red 10 count	12/Ctn
TYC31142386	Devon® 1920 needlecounter bladeguard single foam magnet red 20 count	12/Ctn
TYC31142360	Devon® 1922 needlecounter bladeguard foam adhesive red 20 count	12/Ctn
DRY25-0104	DeRoyal® Needlecounter black double magnet red 30 count	50/Ctn
DRY25-1104Y	DeBlade™ Needlecounter double magnet yellow 30 count	50/Ctn
DRY25-1502	DeBlade™ Needlecounter foam adhesive red 20 count	20/Ctn
DRY25-1401Y	DeBlade™ Needlecounter foam magnet yellow 20 count	20/Ctn
DRY25-1402Y	DeBlade™ Needlecounter foam magnet yellow 10 count	25/Ctn
DRY25-1504	DeBlade™ Needlecounter foam/adhsve 40ct 40 count/80 capacity red	12/Ctn
DRY25-1701	DeBlade™ Needlecounter magnet foam strips red 30 count	50/Ctn
DRY25-0112Y	DeBlade™ Needlecounter double magnet yellow 20 count	65/Ctn
DRY25-0504	DeBlade™ Needlecounter foam adhesive 40ct No, 40 count/80 capacity red	12/Ctn
DRY25-0802	DeBlade™ Needlecounter foam strip red 15 count sterile	50/Ctn
DRY25-1801YNS	DeBlade™ Needlecounter foam strips 30ct Non-Sterile Yellow BNS	100/Ctn
DRY25-1801Y	DeBlade™ Needlecounter foam strips yellow 30 count sterile	50/Ctn
DRY25-0103Y	DeBlade™ Needlecounter single magnet yellow 10 count sterile	50/Ctn
DRY25-0105	DeBlade™ Needlecounter single magnet with clear top red 30 count sterile	27/Ctn



Surgical Consumables continued

Instrument Wipes & Anti-Fog

Item Code	Description	UOM
DRY28-0102	DeRoyal® Defogger™ flip lid w/foam pad Sterile Cylindrical bottle, flip top lid	20/Ctn
DRY28-0101	DeRoyal® Defogger™ w/foam pad Sterile Oval bottle, screw top lid	20/Ctn
SIO400751	ELVIS Dr. Fog® Anti-Fog bottle w/sponge Sterile	10/Ctn
SIO400750	ELVIS Dr. Fog® Anti-Fog pre-treated sponge Sterile	10/Ctn
DRY25-201	DeRoyal® Instrument wipe delicate Sterile scratch-free micro tip wipe	20/Pack
TYC31142535	Devon® Micro tip wipe 3300	25/Ctn



Light Handle Covers

Item Code	Description	UOM
DRY26-011	DeRoyal® Surgiclick™ Light handle cover Sterile single pack	60/Ctn
DRY26-012	DeRoyal® Surgiclick™ Light handle cover Sterile twin pack	300/Ctn
TYC571711	Devon® Lite Glove Non-Sterile	1000/Ctn
TYC31140208	Devon® Lite Glove Sterile	25/Ctn



Patient Preparation

Item Code	Description	UOM
BIP834779	BIC® Facial Shaver Pkt15	10/Ctn
BIP899416	BIC® Medical prep shaver (20pks × 12s)	240/Ctn
BIP7119	BIC® Sensitive skin razors bulk	1000/Ctn
SMSMH0806	KAI Skin Prep shaver - double sided Non-Sterile PTFE coated razor blades single-use shaver	100/Ctn



promed® Sterile & Non-Sterile Holloware

Item Code	Description	UOM
PMD6570102	promed® Gallipot 60mL Sterile	1000/Ctn
PMD5631	promed® Gallipot 150ml Sterile	50/Ctn
PMD5630	promed® Gallipot / Kidney Dish Pack Sterile	50/Ctn
PMD500SJ	promed® Jug 500mL Sterile	12/Ctn
PMD5801	promed® Jug 1.2L Sterile	10/Ctn
PMD6562102	promed® Kidney Dish 700mL Clear Sterile	50/Ctn
PMD6564102	promed® Kidney Dish 700mL Green Sterile	50/Ctn
PMD6563102	promed® Kidney Dish 700mL Yellow Sterile	50/Ctn
PMD1507	promed® Tray Large Blue 4 Compartment Non-Sterile	900/Ctn



Surgical Consumables continued

promed® Sterile Instruments / Instrument Sets

The promed® range of Sterile Single Use instruments are disposable, cost effective, quality surgical instruments. They come sterile and individually packaged in an easy to open peel pouch. To easily identify as 'single use' they feature blue colour coding on handles as well as single use symbol.

Item Code	Description	UOM
PMD155	promed® Adson tissue forceps sterile single use instrument 12.5cm 1:2 teeth	50/Ctn
PMD156	promed® Adson tissue forceps sterile single use instrument 12.5cm plain serrated	50/Ctn
PMD221	promed® Derf needle holders sterile single use instrument 12.5cm	50/Ctn
PMD110	promed® Dissecting scissors sterile single use instrument 12.5cm sh/sh straight	50/Ctn
PMD151	promed® Dressing forceps sterile single use instrument 12.5cm fine	50/Ctn
PMD150	promed® Dressing forceps sterile single use instrument 12.5cm standard	50/Ctn
PMD100	promed® Dressing scissors sterile single use instrument 12.5cm sh/bl straight	50/Ctn
PMD6843	promed® Instrument pack - thyroid tray Sterile	50/Ctn
PMD6613	promed® Instrument Set Sterile	15/Ctn
PMD121	promed® Iris scissors sterile single use instrument 11.5cm curved	50/Ctn
PMD120	promed® Iris scissors sterile single use instrument 11.5cm straight	50/Ctn
PMD181	promed® Magill sponge forceps sterile single use instrument 16cm	50/Ctn
PMD182	promed® Magill sponge forceps sterile single use instrument 20cm	50/Ctn
PMD183	promed® Magill sponge forceps sterile single use instrument 24cm	50/Ctn
PMD131	promed® Mayo Scissors sterile single use instrument 12.5cm curved	50/Ctn
PMD130	promed® Mayo Scissors sterile single use instrument 12.5cm straight	50/Ctn
PMD133	promed® Mayo Scissors sterile single use instrument 14cm curved	50/Ctn
PMD132	promed® Mayo Scissors sterile single use instrument 14cm straight	50/Ctn
PMD220	promed® Mayo-Hegar needle holders sterile single use instrument 15cm	50/Ctn
PMD6844	promed® Micro suture set sterile single use	50/Ctn
PMD6617	promed® Minor Procedure Instruments Pack sterile single use instrument	20/Ctn
PMD4002	promed® Mosquito Forcep 12.5cm	150/Ctn
PMD6812	promed® Suture Tray Sterile	25/Ctn
PMD180	promed® Rampley sponge forceps sterile single use instrument 18cm	50/Ctn
PMD240	promed® Scalpel handle #3 sterile single use instrument	50/Ctn
LVSWHFSW125S	promed® Spencer Wells forceps 12.5cm straight	1/Pack
PMD100	promed® Surgical scissors sterile single use instrument 12.5cm sharp/blunt	50/Ctn
PMD200	promed® Suture Pack 1: scissors, forceps, needle holder	50/Ctn
PMD6745	promed® Suture Pack sterile	10/Ctn
PMD201	promed® Suture Pack 2: scissors, forceps × 2, needle holder	30/Ctn
PMD6772	promed® Suture Pack sterile	20/Ctn
PMD202	promed® Suture Removal pack 1: scissors & forceps	50/Ctn
PMD4202	promed® Suture Set sterile	25/Ctn
PMD3011	promed® Tracheal Dilator sterile 14cm	50/Ctn
PMD230	promed® Tubing clamps sterile single use instrument 16cm	50/Ctn
PMD140	promed® Umbilical cord scissors sterile single use instrument	50/Ctn



Surgical Consumables continued

promed® Pacemaker Instrument Pack

Item Code: PMD6952

UOM: 10/Ctn

Contains:

2 × Prep Swab Stick	1 × Forcep Adson Toothed
2 × Forcep Artery Mosquito Curved 12.5cm	2 × Criles Straight 14cm
1 × Forcep Spencer Wells Curved	1 × Gauze Swabs 10cm × 10cm 8ply (Pkt5)
1 × Ligature Needle	1 × Metzenbaum Scissor Curved 16cm
1 × Needle Holder	2 × Forceps Plastic
1 × Scissor Iris Sharp/Sharp 11.5cm	1 × Towel Clamp
1 × Mollison Retractors Sharp 16cm	1 × Kidney Dish Clear 700ml
2 × Gallipot Clear 150ml	1 × Forcep Dressing
1 × 3-Compartment Tray	1 × Bowl 1L
1 × Wrap Heavy Duty SMS 130cm × 150cm	

promed® Patient and Staff Gowns

Item Code	Description	UOM
PMD182060	promed® Cover Gown, Blue, X-Small	40/Ctn
PMD182070	promed® Cover Gown, Blue, Small	40/Ctn
PMD182080	promed® Cover Gown, Blue, Medium	40/Ctn
PMD182090	promed® Cover Gown, Blue, Large	40/Ctn
PMD182095	promed® Cover Gown, Blue, X-Large	40/Ctn
PMD181060	promed® Cover Gown, White, X-Small	40/Ctn
PMD181070	promed® Cover Gown, White, Small	40/Ctn
PMD181080	promed® Cover Gown, White, Medium	40/Ctn
PMD181090	promed® Cover Gown White Large	40/Ctn
PMD181095	promed® Cover Gown, White, X-Large	40/Ctn
PMD185060	promed® Patient Gown, Dark Blue, Children (XX-Small)	50/Ctn
PMD185070	promed® Patient Gown, Dark Blue, Small	50/Ctn
PMD185080	promed® Patient Gown, Dark Blue, Medium	50/Ctn
PMD185090	promed® Patient Gown, Dark Blue, Large	50/Ctn
PMD185095	promed® Patient Gown, Dark Blue, X-Large	50/Ctn
PMD185096	promed® Patient Gown, Dark Blue, Bariatric (XX-Large)	50/Ctn
PMD180070T	promed® Scrub Suit Dark Blue Top Small	40/Ctn
PMD180080T	promed® Scrub Suit Dark Blue Top Medium	40/Ctn
PMD180090T	promed® Scrub Suit Dark Blue Top Large	40/Ctn
PMD180095T	promed® Scrub Suit Dark Blue Top X-Large	40/Ctn
PMD180096T	promed® Scrub Suit Dark Blue Top XX-Large	40/Ctn
PMD180070B	promed® Scrub Suit Pants Dark Blue Small	40/Ctn
PMD180080B	promed® Scrub Suit Pants Dark Blue Medium	40/Ctn
PMD180090B	promed® Scrub Suit Pants Dark Blue Large	40/Ctn
PMD180095B	promed® Scrub Suit Pants Dark Blue X-Large	40/Ctn
PMD180096B	promed® Scrub Suit Pants Dark Blue XX-Large	40/Ctn
PMD187060	promed® Warming Jacket, Blue, X-Small	40/Ctn
PMD187070	promed® Warming Jacket, Blue, Small	40/Ctn
PMD187080	promed® Warming Jacket, Blue, Medium	40/Ctn
PMD187090	promed® Warming Jacket, Blue, Large	40/Ctn
PMD187095	promed® Warming Jacket, Blue, X-Large	40/Ctn



Ultrasound Devices

Thermasonic® Gel Warmer

The NEW single bottle and two re-designed multi-bottle gel warmers rapidly warm ultrasound gel and consistently maintain temperature. Durable construction, suitable for wall mounting or counter top use and each unit comes with a two-year warranty. Options available include: Single Bottle with indicator, 3-Bottle holder with LED and 3-Bottle holder with LCD.

Features:

- Powerful heating element uniformly heats bottles of gel
- High-density, impact-resistant composite
- Microprocessor control maintains constant temperature
- Adjustable temperature settings for multi-bottle units
- “Thermal cutoff” provides protection from overheating
- Illuminated power ON/OFF switch on multi-bottle units
- Light-weight and portable design
- Heat Indicator Lamp on single-bottle unit
- Low electrical demand

Item Code	Description	UOM
PKR82-01-20CE	Thermasonic® Single-bottle gel warmer 230V	EA
PKR82-03-20CE	Thermasonic® Multi-bottle gel warmer 230V with LEDs	EA
PKR83-03-20LCD	Thermasonic® Multi-bottle gel warmer 230Vwith LCD display	EA



Ultrasound Devices continued

signosRT® - Portable Ultrasound Device

POCKET-SIZED AND PERSONAL – INSTANT INSIGHT at the POINT-OF-CARE

signosRT®, the world's smallest and most affordable completely portable ultrasound system, now makes ultrasound technology easy to use – for all users from novice to experienced. It is small enough to carry in your pocket or wear like a stethoscope.

Providing outstanding image quality, this unique palm-sized device weighs less than 400 grams (14 ounces) and is ideal for quick and simple point-of-care applications.

Some popular applications of the signosRT® system include

- eFAST scan
- Lung imaging - pleural fluid, pneumothoraces
- Post-void residual bladder volume (automated measurements)
- Basic pregnancy assessments – confirmation scans, FHR, fetal position
- Basic internal organ imaging – gross gall bladder pathology.
- Abdominal aortic aneurysm (AAA) Screening
- Localisation of free fluid in the abdomen and lungs
- Core muscle assessment – TA, pelvic floor
- Basic veterinary applications

Item Code	Description	UOM
SIGP03777	signosRT® Ultrasound Package S3-5 Transducer, incl. Software & Cables	1/Pack
SIGP03611	signosRT® S3-5 Transducer	1/Pack
SIGRT2YWTY	signosRT® 2 Year Extended Warranty	-
SIGP03620	signosRT® Data connection cable for PC, Laptop, Tablet	1/Pack

Easy access to USB port and headphone jack, touch screen stylus and 4GB microSD data card (included) which enables you to store 10,000+ images.



Convenient USB connection enables removable probe to be directly connected to your Windows based PC, laptop or tablet running the accompanying SigViewer software.

High resolution touch screen display shows accurate images and enables you to easily and quickly access a wide range of applications – Presets, Scanning Modes, Depth, Gain, Measurements, Patient Information, Text Annotation and Voice Recording – for fast clinical review.



The 'S3-5' is suitable for all applications of the signosRT®, providing higher resolution images from superficial areas through to 18cm.

Ultrasound Consumables

Aquasonic 100® Ultrasound Transmission Gel

World standard for medical ultrasound transmission, Aquasonic 100 is the widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required.

Features:

- Correct acoustics for the broad range of frequencies used
- Completely aqueous, will not stain clothing or damage transducers
- Unique “can’t be copied” formula is hypoallergenic, bacteriostatic, non-sensitising and non irritating
- Used and recommended by leading manufacturers of medical ultrasound equipment worldwide
- No formaldehyde
- Not a spermicide

Item Code	Description	UOM
PKR01-20	Aquasonic 100® Ultrasound Gel 20g Non-Sterile Single Use	100/Pack
PKR01-02	Aquasonic 100® Ultrasound Gel 60g	12/Pack
PKR01-08	Aquasonic 100® Ultrasound Gel 250ml	12/Pack
PKR01-50	Aquasonic 100® Ultrasound Gel 5L	1/Pack



Sterile Aquasonic 100® Ultrasound Transmission Gel

World Standard for Sterile ultrasound gel. Sterile Aquasonic® 100 is the recommended for diagnostic and therapeutic ultrasound applications when sterility is indicated.

Features:

- Overwrapped, sterilised foil pouches
- Acoustically correct for the broad range of frequencies used in medical ultrasound
- Non-irritating, hypoallergenic
- Completely aqueous, will not stain clothing or damage equipment
- Not a spermicide

Item Code	Description	UOM
PKR01-01	Aquasonic 100® Ultrasound Gel 20g single use sterile pack	48/Box



Ultrasound Consumables continued

Aquasonic Clear® Ultrasound Gel

Clear ultrasound gel. Based on the same “can’t be copied” formula of the World Standard, Aquasonic 100 Ultrasound Transmission Gel, without fragrance or colour. For diagnostic and therapeutic medical ultrasound examinations.

Features:

- No fragrance or colour
- Correct acoustics for the broad range of frequencies used
- Hypoallergenic, bacteriostatic, non-irritating
- Water soluble
- No formaldehyde
- Not a spermicide

Item Code	Description	UOM
PKR03-08	Aquasonic Clear® Ultrasound Gel 250mL	12/Pack
PKR03-50	Aquasonic Clear® Ultrasound Gel 5L	1/Box



Aquaflex® Ultrasound Gel Pad

A 2cm × 9cm aqueous, flexible, disposable ultrasound standoff for use in difficult to visualise and near field areas, or when gels alone won't do. Disposability eliminates the risk of cross contamination in diagnostic and therapeutic ultrasound procedures.

Features:

- Requires no gel, except in prolonged procedures
- Requires only gentle, even pressure
- Hypoallergenic, Bacteriostatic
- Cushions sensitive areas
- Easy to use
- Disposable

Item Code	Description	UOM
PKR04-02	Aquaflex® Ultrasound Gel pad 2 × 9cm	6/Pack



Transeptic® Cleaning Solution

Multi-purpose spray for professional use when cleansing ultrasound transducer/probe surfaces between procedures. Does not cause transducer membrane swelling.

Features:

- Convenient clear spray dispenser
- Does not contain phenol or glutaraldehyde
- Accepted for use by manufacturers of ultrasound equipment

Item Code	Description	UOM
PKR09-25	Transeptic® Cleansing solution 250ml	12/Pack



Ultrasound Consumables continued

Scan® Ultrasound Gel

Our exclusive dual-viscosity formula allows SCAN to be applied as a gel and scanned as a fluid. For medical ultrasound procedures when a less viscous gel is desired.

Features:

- Accoustically correct, provides superior transmission of ultrasound waves
- Completely aqueous, will not stain clothing, damage equipment or irritate skin
- Hypoallergenic, bacteriostatic with no formaldehyde

Item Code	Description	UOM
PKR11-28S	Scan® Ultrasound gel; 4 × 3.8L of SCAN, 2 refillable dispensers and 1 pump	4/Pack



Polysonic® Ultrasound Lotion

Multi-purpose ultrasound lotion with superior coupling efficiency. Recommended for all medical ultrasound procedures where a lotion is preferred.

Features:

- Rich moisturising formula will not dry or irritate skin
- Accoustically correct for the broad range of frequencies used in medical ultrasound
- Will not damage equipment or clothing
- Comfortable & pleasing to patient and ultrasound practitioner
- Hypo-allergenic, bacteriostatic
- No salt, no formaldehyde

Item Code	Description	UOM
PKR21-08	Polysonic® Ultrasound lotion 250mL	12/Pack
PKR21-28	Polysonic® Ultrasound lotion 3.8L	4/Pack



Redux® Products

The Redux® gel is non-abrasive, water-based ultrasound consumable for routine ECG and difibrillation.

- Highly conductive, water soluble, non-staining, non-gritty

Redux® cream is mainly used for routine ECG and not recommended for defibrillation

- Cosmetic quality, pleasing to patients
- Non-irritating
- Highly conductive electrode cream

The Redux® Paste is a water-based, abrasive gel for short term monitoring, ECG and defibrillation; highly-conductive and water soluble.

Item Code	Description	UOM
PKR66-04	Redux® Electrolyte cream 114g 4oz tube	10/Pack
PKR65-04	Redux® Electrolyte gel 114g 4oz tube	10/Pack
PKR67-05	Redux® Electrolyte paste 142g 5oz tube	10/Pack



Ultrasound Consumables continued

Signacreme®

Highly conductive cream electrolyte for use in electromedical procedures where a cosmetic quality cream is preferred.

- Cosmetic quality, pleasing to patients
- Non-irritating, hypoallergenic, bacteriostatic
- Not recommended for defibrillation
- Conveniently packaged in refillable see-thru dispenser with flip-top cap

Item Code	Description	UOM
PKR17-05	Signacreme® Electrode cream 142g	12/Pack



Signagel® Electrode Gel

Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback and EMG.

- Highly conductive
- Water soluble, non-staining, non-gritty
- Hypoallergenic, bacteriostatic
- See-thru tube with stand-up cap lets you know how much Signagel® is left

Item Code	Description	UOM
PKR15-60	Signagel® Electrode gel 60g	12/Pack
PKR15-25	Signagel® Electrode gel clear tube 250g	12/Pack



Signaspray® Electrode Solution And Skin Preparation

Highly conductive, economical spray electrolyte and non-gritty skin prep.

- The ideal electrolyte for use with muscle stimulation equipment
- Bacteriostatic, hypoallergenic, non-staining
- Clear bottle lets you know how much Signaspray® is left
- Superior non-gritty skin prep
- No residue build-up
- Easy to apply, easy to clean up
- Not recommended for defibrillation

Item Code	Description	UOM
PKR18-25	Signaspray® Electrode solution-skin prep 250mL	12/Pack



Probe Covers

Eclipse® / Eclipse® 3D - Probe Cover - Latex Free

A soft, latex-free, disposable probe cover; one that is designed for standard ultrasound probes/transducers, and the other is designed specifically for 3D/4D ultrasound probes/transducers. Both are pregelled inside with Aquasonic 100 Ultrasound transmission Gel.

Features:

- Pre-gelled inside with Aquasonic 100 Ultrasound Transmission Gel
- Latex-free, non-irritating, non-sensitising
- Only one probe cover is required for each procedure

Item Code	Description	UOM
PKR38-03	Eclipse® 3D Probe Cover latex free; 83 × 43 × 241mm	100/Pack
PKR38-01	Eclipse® Probe Cover latex free 64 × 44 × 241mm	100/Pack



Civco Probe Covers

Item Code	Description	UOM
PHL610-249	Civco Non-Sterile Rolled Latex cover 3.5 × 20cm	100/Pack
PHL610-362	Civco Non-Sterile Civ-flex cover 14 × 61cm (3D)	50/Pack
PHL610-208	Civco Sterile polyethylene cover 15.2 × 83.8cm	24/Pack
PMD5737	Civco Sterile Probe Cover Kit – includes sterile Aquasonic Gel Sachet & two rubber bands	25/Pack

promed® Sterile Ultrasound Probe Covers

Item Code	Description	UOM
PMD5370	promed® Probe cover with rubber band Sterile 12cm × 100cm	50/Ctn
PMD8225102	promed® Probe cover Sterile 10cm × 132cm	50/Ctn
PMD8224102	promed® Probe cover Sterile 10cm × 100cm	50/Ctn
PMD5336.0	promed® Probe cover Sterile 30cm × 90cm	100/Ctn
PMD5337.0	promed® Probe cover Sterile 44cm × 60cm	100/Ctn
PMD5337	promed® Probe cover Sterile 45cm × 60cm	100/Ctn
PMDA113	promed® Ultrasound Probe Cover Sterile 15cm × 225cm	75/Ctn
PMDU2115	promed® Ultrasound Probe Cover Sterile 20cm × 120cm	100/Ctn
PMDU6103	promed® Ultrasound Probe Cover Sterile 20.8cm × 50cm	100/Ctn
PMDU6105	promed® Ultrasound Probe Cover Sterile 22.8cm × 100cm	100/Ctn





X-Ray Consumables

Includes Image Printing, X-Ray Envelopes, Labels and Stickers:

We stock many radiology staples such as film stickers, film bag labels and X-Ray Markers (metal stencils).

We also offer made to order products so that you may customise your labels or stencils with radiographer initials or specific wording.

Image Printing

Item Code	Description	UOM
JCMMUPP110HD	Print Paper HD No2 Sony	5/Pack
JCMMUPP110HG	Print Paper Roll High Gloss Sony	5/Pack
JCMMUPP210HD	Print Paper Roll for Sony Model UP-980CU	5/Pack
JCMMUPC-510	Sony Colour Printing Pack Ensemble	200/Pack
JCMMUPC-21L	Sony Print Pack Film	4/Pack
JCMMK61S	Mitsubishi thermal paper	4/Pack

X-Ray Envelopes

Item Code	Description	UOM
JCMS01-140267	Gold Envelopes 444 × 368mm	250/Pack
JCMS08-140216	Grey envelopes 317 × 266mm	250/Pack
JCMS07-140219	Grey envelopes 444 × 368mm	250/Pack
JCMSR706-140193	OPG Gold envelopes 355 × 152mm	250/Pack
JCMS15-140217	White envelopes 317 × 266mm	250/Pack
JCMS14-140220	White envelopes 444 × 368mm	250/Pack

X-Ray Consumables continued

X-Ray Markers (Metal stencils)

Item Code	Description	UOM
JCBS-01	Stencil 1 letter/digit 10×10mm	1/Pack
JCBS-02	Stencil 2 letter/digit 15×10mm	1/Pack
JCBS-03	Stencil 3 letter 25×10mm	1/Pack
JCBS-04	Stencil 4 letter 30×10mm	1/Pack
JCBS-05	Stencil 5 letter 40×10mm	1/Pack
JCBS-06	Stencil 6 letter 45×10mm	1/Pack
JCBS-07	Stencil 7 letter 50×10mm	1/Pack
JCBS-08	Stencil 8 letter 60×10mm	1/Pack
JCBS-09	Stencil 9 letter 65×10mm	1/Pack
JCBS-10	Stencil 10 letter 70×10mm	1/Pack
JCBS-11	Stencil 11 letter 80×10mm	1/Pack
JCBS-12	Stencil 12 letter 85×10mm	1/Pack
JCBS-13	Stencil 13 letter 90×10mm	1/Pack
JCBS-14	Stencil 14 letter 100×10mm	1/Pack
JCBS-15	Stencil 15 letter 105×10mm	1/Pack
JCBS-16	Stencil 16 letter 110×10mm	1/Pack
JCBS-17	Stencil 17 letter 120×10mm	1/Pack
JCBS-18	Stencil 18 Letter 125×10mm	1/Pack
JCBS-19	Stencil 19 letter 130×10mm	1/Pack
JCBS-20	Stencil 20 letter 140×10mm	1/Pack
JCBS-23	Stencil 23 letter 160×10mm	1/Pack
JCBS-26	Stencil 26 letter 180×10mm	1/Pack
JCBS-1I	Stencil R or L with 1 letter/digit 20×20mm	1/Pack
JCBS-2I	Stencil R or L 2 initials 20×20mm	1/Pack
JCBS-3I	Stencil R or L with 3 vertical initials 25×25mm	1/Pack
JCBS-4I	Stencil R or L with 4 vertical initials 30×30mm	1/Pack
JCBS-5I	Stencil R/L with 5 initials	1/Pack
JCBS-6I	Stencil R/L with 6 initials	1/Pack
JCBS-CC	Stencil CC 15×10mm	1/Pack
JCBS-POS-MICT	Stencil POST MICT 60×10mm	1/Pack
JCBS-EXTROT	Stencil EXT ROT 50×10mm	1/Pack
JCBS-PRELIM	Stencil PRELIM 40×10mm	1/Pack
JCBS-DECUBITUS	Stencil DECUBITUS 60×10mm	1/Pack
JCBS-SUPINE	Stencil SUPINE 40×10mm	1/Pack
JCBS-MOBILE	Stencil MOBILE 40×10mm	1/Pack
JCBS-CLOSED	Stencil CLOSED 45×10mm	1/Pack
JCBS-ROUNDL	Stencil round brass L	1/Pack
JCBS-ROUND R	Stencil round brass R	1/Pack
JCBS-01L	Stencil LEFT 10×10mm	1/Pack
JCBS-01R	Stencil RIGHT 10×10mm	1/Pack
JCBS-1NL	Stencil L with 1 digit 20×20mm	1/Pack
JCBS-1NR	Stencil R with 1 digit 20×20mm	1/Pack
JCBS-2NL	Stencil L with 2 digits 20×20mm	1/Pack
JCBS-2NR	Stencil R with 2 digits 20×20mm	1/Pack

X-Ray Consumables continued

Item Code	Description	UOM
JCBS-3IH	Stencil R or L with 3 horizontal initials 25×25mm	1/Pack
JCBS-3NL	Stencil L with 3 vertical digits 25×25mm	1/Pack
JCBS-3NR	Stencil R with 3 vertical digits 25×25mm	1/Pack
JCBS-EXP	Stencil EXP 25×10mm	1/Pack
JCBS-EXT	Stencil EXT 25×10mm	1/Pack
JCBS-MAG	Stencil MAG 25×10mm	1/Pack
JCBS-MLO	Stencil MLO 25×10mm	1/Pack
JCBS-OBL	Stencil OBL 25×10mm	1/Pack
JCBS-SET	Stencil set engraved 0-9 10×10mm	1/Pack
JCBS-COMPRESSION	Stencil COMPRESSION 75×10mm	1/Pack
JCBS-MAMMOSET	Stencil set MLO CC LAT CLEO MAG - 1 of each	1/Pack
JCBS-OBLIQUE	Stencil OBLIQUE 45×10mm	1/Pack
JCBS-01L-20	Stencil L 20×20mm	1/Pack
JCBS-01R-20	Stencil R 20×20mm	1/Pack
JCBS-LLAT-OBL	Stencil LLAT-OBL 60×10mm	1/Pack
JCBS-CONTROL	Stencil CONTROL 50×10mm	1/Pack
JCBS-LATERAL	Stencil LATERAL 50×10mm	1/Pack
JCBS-AP SITTING	Stencil AP SITTING	1/Pack
JCBS-CLEO	Stencil CLEO 30×10mm	1/Pack
JCBS-INSP	Stencil INSP 25×10mm	1/Pack
JCBS-LEFT	Stencil LEFT 30×10mm	1/Pack
JCBS-OPEN	Stencil OPEN 30×10mm	1/Pack
JCBS-PREM	Stencil PREM 30×10mm	1/Pack
JCBS-PRONE	Stencil PRONE 35×10mm	1/Pack
JCBS-EXPIRATION	Stencil EXPIRATION 65×10mm	1/Pack
JCBS-RIGHT	Stencil RIGHT 35×10mm	1/Pack
JCBS-RLAT-OBL	Stencil RLAT-OBL 60×10mm	1/Pack
JCBS-2I-10MM	Stencil R or L 2 initials 10×10mm	1/Pack
JCBS-ALPHA	Stencils set brass A-Z 10×10mm each	1/Pack
JCBS-SET20	Stencil set engraved 0-9 20×20mm	1/Pack
JCBS-EXTERNALROT	Stencil EXTERNAL ROTATION 110×10mm	1/Pack
JCBS-RELEASE	Stencil RELEASE 50×10mm	1/Pack
JCBS-ARROW	Stencil pointed arrow 35×10mm	1/Pack

X-Ray Consumables continued

X-Ray Stickers and Labels

Item Code	Description	UOM
JCBLBOY	I am a good boy blue 300/Pack	300/Pack
JCBLGIRL	I am a good girl pink 300/Pack	300/Pack
JCBLDOLLYRAY	I sat still x-ray doll green/white	300/Pack
JCBLPUPXRAY	I sat still x-ray doll red/black	300/Pack
JCBLSMILES	Smiley sun/flowers orange/yellow	400/Pack
JCBL4238	Allergy 30×11mm black/fluoro red	500/Pack
JCBL4223	Caution patient same name black/fluoro orange	250/Pack
JCBL4335	Confidential refer doctor 120×35mm	250/Pack
JCBL4217	Confidential refer doctor 63×20mm	250/Pack
JCBL4229	CT scans performed 48×21mm blue/white	250/Pack
JCBL4130	Duplicate black/white	500/Pack
JCBL5009	Echo cardiography performed black/fluoro red	200/Pack
JCBL4184	Erect red/white	500/Pack
JCBL4180	Inpatient red/white	500/Pack
JCBL7515	Label 25×13mm white-plain	250/Pack
JCBL4174	Left black/white	500/Pack
JCBL4279	Medical imaging 45×11 white/blue	250/Pack
JCBL4185	Mobile red/white	500/Pack
JCBL4150	Oblique black/white	500/Pack
JCBL4216	Patient films 68×30mm black/fluoro red	200/Pack
JCBL4224	Please do not destroy black/red	250/Pack
JCBL4281	Previous films 68×30mm black/fluoro red	250/Pack
JCBL4172	Square L black/white	1000/Pack
JCBL4173	Square R black/white	1000/Pack
JCBL4182	Supine Red/White	500/Pack
JCBL4201	U/S label fold over black/pale blue	1000/Pack
JCBL4250	Ultrasound 48×21mm black/fluoro yellow	250/Pack
JCBL4230	Ultrasound performed 48×21mm blue/white	250/Pack
JCBL4221	Unreported pls rtn X-Ray 48×21mm black/red	250/Pack
JCBL4221RE	Unreported pls rtn X-Ray 48×21mm removable	250/Pack
JCBL4236	Urgent 62×20mm white/red	250/Pack
JCBL4151	With contrast black/white	500/Pack
JCBL4215	X-Ray films do not bend 48×21mm	250/Pack
JCBL4369	X-Rays taken 48×21mm blue/white	250/Pack
JCBL4301	Mammography 48×21mm black/fluoro red	250/Pack
JCBL4227	Mobile unit/erect 67×31mm black/white	250/Pack
JCBL6000-12	12/12 year label black/blue	250/Pack
JCBL6000-13	12/13 year label fluoro orange & black	250/Pack
JCBL6000-14	13/14 year label black/purple	250/Pack



X-Ray Consumables continued

Item Code	Description	UOM
JCBL00A	Letter A 13mm black/white	1000/Pack
JCBL00B	Letter B 13mm black/white	1000/Pack
JCBL00C	Letter C 13mm black/white	1000/Pack
JCBL00D	Letter D 13mm black/white	1000/Pack
JCBL00E	Letter E 13mm black/white	1000/Pack
JCBL00F	Letter F 13mm black/white	1000/Pack
JCBL00G	Letter G 13mm black/white	1000/Pack
JCBL00H	Letter H 13mm black/white	1000/Pack
JCBL00I	Letter I 13mm black/white	1000/Pack
JCBL00J	Letter J 13mm black/white	1000/Pack
JCBL00K	Letter K 13mm black/white	1000/Pack
JCBL00L	Letter L 13mm black/white	1000/Pack
JCBL00M	Letter M 13mm black/white	1000/Pack
JCBL00O	Letter O 13mm black/white	1000/Pack
JCBL00P	Letter P 13mm black/white	1000/Pack
JCBL00Q	Letter Q 13mm black/white	1000/Pack
JCBL00R	Letter R 13mm black/white	1000/Pack
JCBL00S	Letter S 13mm black/white	1000/Pack
JCBL00T	Letter T 13mm black/white	1000/Pack
JCBL0000	Figure 0 13mm black/white	1000/Pack
JCBL0001	Figure 1 13mm black/white	1000/Pack
JCBL0002	Figure 2 13mm black/white	1000/Pack
JCBL0003	Figure 3 13mm black/white	1000/Pack
JCBL0004	Figure 4 13mm black/white	1000/Pack
JCBL0005	Figure 5 13mm black/white	1000/Pack
JCBL0006	Figure 6 13mm black/white	1000/Pack
JCBL0007	Figure 7 13mm black/white	1000/Pack
JCBL0008	Figure 8 13mm black/white	1000/Pack
JCBL0009	Figure 9 13mm black/white	1000/Pack
JCBL0010	Figure 10 13mm black/white	1000/Pack
JCBL0011	Figure 11 13mm black/white	1000/Pack
JCBL0012	Figure 12 13mm black/white	1000/Pack
JCBL0013	Figure 13 13mm black/white	1000/Pack
JCBL0014	Figure 14 13mm black/white	1000/Pack
JCBL0015	Figure 15 13mm black/white	1000/Pack
JCBL0016	Figure 16 13mm black/white	1000/Pack
JCBL0017	Figure 17 13mm black/white	1000/Pack
JCBL0018	Figure 18 13mm black/white	1000/Pack
JCBL0019	Figure 19 13mm black/white	1000/Pack
JCBL0020	Figure 20 13mm black/white	1000/Pack

X-Ray Consumables continued

Item Code	Description	UOM
JCBL1001	Number 1 - series 1 grey colour	1000/Pack
JCBL1002	Number 2 - series 1 red colour	1000/Pack
JCBL1003	Number 3 - series 1 lemon colour	1000/Pack
JCBL1004	Number 4 - series 1 light green	1000/Pack
JCBL1005	Number 5 - series 1 dark blue	1000/Pack
JCBL1006	Number 6 - series 1 black/white	1000/Pack
JCBL1007	Number 7 - series 1 mauve colour	1000/Pack
JCBL1008	Number 8 - series 1 yellow colour	1000/Pack
JCBL1009	Number 9 - series 1 teal colour	1000/Pack
JCBL1010	Number 10 - series 1 pink colour	1000/Pack
JCBL1011	Number 11 - series 1 grey stripe	1000/Pack
JCBL1012	Number 12 - series 1 red stripe	1000/Pack
JCBL1013	Number 13 - series 1 lemon stripe	1000/Pack
JCBL1014	Number 14 - series 1 green stripe	1000/Pack
JCBL1015	Number 15 - series 1 blue stripe	1000/Pack
JCBL1016	Number 16 - series 1 black on white	1000/Pack
JCBL1017	Number 17 - series 1 mauve stripe	1000/Pack
JCBL1018	Number 18 - series 1 yellow stripe	1000/Pack
JCBL1019	Number 19 - series 1 teal stripe	1000/Pack
JCBL1020	Number 20 - series 1 pink stripe	1000/Pack

3
.....
£

Series 1
JCBL1001-1020


JCBL2001	Number 1 - series 2 grey colour	1000/Pack
JCBL2002	Number 2 - series 2 red colour	1000/Pack
JCBL2003	Number 3 - series 2 lemon colour	1000/Pack
JCBL2004	Number 4 - series 2 light green	1000/Pack
JCBL2005	Number 5 - series 2 dark blue	1000/Pack
JCBL2006	Number 6 - series 2 black/white	1000/Pack
JCBL2007	Number 7 - series 2 mauve colour	1000/Pack
JCBL2008	Number 8 - series 2 yellow colour	1000/Pack
JCBL2009	Number 9 - series 2 teal colour	1000/Pack
JCBL2010	Number 10 - series 2 pink colour	1000/Pack
JCBL2011	Number 11 - series 2 grey stripe	1000/Pack
JCBL2012	Number 12 - series 2 red stripe	1000/Pack
JCBL2013	Number 13 - series 2 lemon stripe	1000/Pack
JCBL2014	Number 14 - series 2 green stripe	1000/Pack
JCBL2015	Number 15 - series 2 blue stripe	1000/Pack
JCBL2016	Number 16 - series 2 black on white	1000/Pack
JCBL2017	Number 17 - series 2 mauve stripe	1000/Pack
JCBL2018	Number 18 - series 2 yellow stripe	1000/Pack
JCBL2019	Number 19 - series 2 teal stripe	1000/Pack
JCBL2020	Number 20 - series 2 pink stripe	1000/Pack

3
.....
3

Series 2
JCBL2001-2020

X-Ray Consumables continued

Item Code	Description	UOM
JCBL3001	Number 1 - series 3 grey colour	1000/Pack
JCBL3002	Number 2 - series 3 red colour	1000/Pack
JCBL3003	Number 3 - series 3 lemon colour	1000/Pack
JCBL3004	Number 4 - series 3 light green	1000/Pack
JCBL3005	Number 5 - series 3 dark blue	1000/Pack
JCBL3006	Number 6 - series 3 black/white	1000/Pack
JCBL3007	Number 7 - series 3 mauve colour	1000/Pack
JCBL3008	Number 8 - series 3 yellow colour	1000/Pack
JCBL3009	Number 9 - series 3 teal colour	1000/Pack
JCBL3010	Number 10 - series 3 pink colour	1000/Pack
JCBL3011	Number 11 - series 3 grey stripe	1000/Pack
JCBL3012	Number 12 - series 3 red stripe	1000/Pack
JCBL3013	Number 13 - series 3 lemon stripe	1000/Pack
JCBL3014	Number 14 - series 3 green stripe	1000/Pack
JCBL3015	Number 15 - series 3 blue stripe	1000/Pack
JCBL3016	Number 16 - series 3 black on white	1000/Pack
JCBL3017	Number 17 - series 3 mauve stripe	1000/Pack
JCBL3018	Number 18 - series 3 yellow stripe	1000/Pack
JCBL3019	Number 19 - series 3 teal stripe	1000/Pack
JCBL3020	Number 20 - series 3 pink stripe	1000/Pack
JCBL3021	Number 21 - series 3 grey colour	1000/Pack
JCBL3022	Number 22 - series 3 red colour	1000/Pack
JCBL3023	Number 23 - series 3 lemon colour	1000/Pack
JCBL3024	Number 24 - series 3 light green	1000/Pack
JCBL3025	Number 25 - series 3 dark blue	1000/Pack
JCBL3026	Number 26 - series 3 black/white	1000/Pack
JCBL3027	Number 27 - series 3 mauve colour	1000/Pack
JCBL3028	Number 28 - series 3 yellow colour	1000/Pack
JCBL3029	Number 29 - series 3 teal colour	1000/Pack
JCBL4412	Page 1 black/white	1000/Pack
JCBL4413	Page 2 black/white	1000/Pack
JCBL4414	Page 3 black/white	1000/Pack
JCBL4415	Page 4 black/white	1000/Pack
JCBL4416	Page 5 black/white	1000/Pack
JCBL4417	Page 6 black/white	1000/Pack
JCBL4418	Page 7 black/white	1000/Pack
JCBL4419	Page 8 black/white	1000/Pack
JCBL4420	Page 9 black/white	1000/Pack
JCBL4421	Page 10 black/white	1000/Pack



Series 3
JCBL3001-3020



Thermo Fisher Scientific provide a range of radiology and angiography consumables and devices, specifically focused on providing positive patient outcomes. We complement our manufactured products with key distributed brands. Thermo Fisher Scientific are your reliable and trusted supply partner to deliver patient care.

ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-292
Email AUInfo@thermofisher.com
Visit us online: www.thermofisher.com.au

©2014 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917