



March into Autumn

with great deals

ThermoFisher
SCIENTIFIC

Thermo Scientific™ Nunc™ cell culture consumables

- 10⁻⁶ sterility assurance level (SAL) as per ISO 11137
- Each batch is tested with 4 cell types
 - 2 cell lines and 1 primary cell culture for monolayer formation
 - 1 cell line for cloning efficiency
- Certified non-pyrogenic
- EasYFlask™ has ergonomically designed 1/3 turn cap
- Multiwell plates have raised well rims to lower risk of cross contamination



Quality from start to finish

Reliability • Reproducibility • Scalability • Traceability • Security • Simplicity

| Item Code | Description | Unit | 1 case | 3+ cases |
|---|--|-------------|----------|----------|
| Cell culture multiwell plates - Nunclon Delta™ cell culture surface - individually wrapped sterile | | | | |
| NUN140675 | Multidish 6 well flat bottom, clear with lid | Case of 75 | \$98.70 | \$85.55 |
| NUN150628 | Multidish 12 well flat bottom, clear with lid | Case of 75 | \$109.95 | \$95.35 |
| NUN142475 | Multidish 24 well flat bottom, clear with lid | Case of 75 | \$91.90 | \$79.65 |
| NUN150687 | Multidish 48 well flat bottom, clear with lid | Case of 75 | \$120.95 | \$104.85 |
| NUN167008 | MicroWell™ 96 well flat bottom, clear with lid | Case of 50 | \$58.70 | \$50.90 |
| Cell culture flasks - Nunclon Delta cell culture surface - sterile | | | | |
| NUN156367 | 25cm ² EasYFlask with filter cap | Case of 200 | \$186.40 | \$161.55 |
| NUN156499 | 75cm ² EasYFlask with filter cap | Case of 100 | \$146.55 | \$127.10 |
| NUN159910 | 175cm ² EasYFlask with filter cap | Case of 30 | \$78.15 | \$67.75 |

Finnpipette™ F1 good laboratory pipetting (GLP) kits

Everything you want in your everyday, go-to pipette

- Silent, effortless volume adjustment
- Lightweight design with lowered pipetting forces reduces risk of RSI
- 5 year warranty



| Item Code | Description | Unit | Price |
|-------------|--|------|----------|
| FNP4700850N | Kit 1: 1-10 µL micro, 10-100 µL and 100-1000 µL pipettes | EA | \$750.00 |
| FNP4700860N | Kit 2: 0.2-2 µL micro, 2-20 µL, 20-200 µL and 100-1000 µL pipettes | EA | \$870.00 |
| FNP4700865N | Kit 3: 10-100 µL, 100-1000 µL and 1-10 mL pipettes | EA | \$750.00 |
| FNP4701060N | Kit 4: 2-20 µL, 20-200 µL and 100-1000 µL pipettes | EA | \$750.00 |

15% off list price on all stepper tips

Finnpipette Stepper Pipette

- Easy to use stepper pipette up to 45 dispensing
- Light weight ergonomic design
- Simple easy to set volume adjustment
- Volume range 10-5000 µL

Stepper Pipette

| Item Code | Description | Unit | Price |
|-----------------------------------|---|--------|----------|
| FNP4540000 | Finnpipette stepper 10-5000 µL | EA | \$420.00 |
| Stepper tips - sterile | | | |
| FNP9404173 | Finntip for Finnpette stepper 0.5 mL sterile | 50pcs | |
| FNP9404183 | Finntip for Finnpette stepper 1.25 mL sterile | 50pcs | |
| FNP9404193 | Finntip for Finnpette stepper 2.5 mL sterile | 50pcs | |
| FNP9404203 | Finntip for Finnpette stepper 5 mL sterile | 25pcs | |
| FNP9404213 | Finntip for Finnpette stepper 12.5 mL sterile | 25pcs | |
| FNP9404223 | Finntip for Finnpette stepper 25 mL sterile | 10pcs | |
| FNP9404233 | Finntip for Finnpette stepper 50 mL sterile | 10pcs | |
| Stepper tips - non-sterile | | | |
| FNP9404170 | Finntip for Finnpette stepper 0.5 mL | 100pcs | |
| FNP9404180 | Finntip for Finnpette stepper 1.25 mL | 100pcs | |
| FNP9404190 | Finntip for Finnpette stepper 2.5 mL | 100pcs | |
| FNP9404200 | Finntip for Finnpette stepper 5 mL | 50pcs | |
| FNP9404210 | Finntip for Finnpette stepper 12.5 mL | 50pcs | |
| FNP9404220 | Finntip for Finnpette stepper 25 mL | 20pcs | |
| FNP9404230 | Finntip for Finnpette stepper 50 mL | 10pcs | |

*Contact customer service on 1300-735-292 for Finntip prices



Thermo Scientific™ Heratherm™ compact microbiological incubators

Ideal for a personalized workspace or crowded labs with small sample volumes

- Temperatures at or below ambient: 17° to 40°C with high temperature accuracy
- Save space with extremely small footprint
- 2 year warranty



| Item Code | Description | Unit | Price |
|-------------|-----------------------------------|------|------------|
| HEA50125882 | IMC18 Heratherm compact incubator | EA | \$1,190.00 |

Dependable reproducible results, flexible, easy use

Thermo Scientific mixing and temperature control drybaths

Maximise lab space and versatility with our compact programmable drybaths

- Choice of heating/mixing, heating/cooling or heating/mixing cooling
- Excellent temperature uniformity across the block
- Easily interchangeable blocks for tubes and plates
- Digital control for easy setup and current status of run
- 2 year warranty



| Item Code | Description | Unit | Price |
|-------------|--|------|------------|
| BAR88880030 | Heating cooling drybath - blocks additional | EA | \$2,040.15 |
| BAR88880028 | Heating shaking drybath - blocks additional | EA | \$2,098.10 |
| GRA13687717 | Heating shaking cooling dry bath with 1.5 ml block | EA | \$3,247.50 |
| GRA13687721 | Heating shaking cooling dry bath with 96-well microplate block | EA | \$3,247.50 |

For increased productivity with reduced muscle effort

Ansell Microflex® 93-833 certified ergonomic gloves

- Certified as ergonomically designed to reduce muscle effort and hand fatigue
- More comfortable and lighter weight than most other examination gloves
- Up to 60% stronger than comparable leading nitrile gloves, providing reliable durability



| Item Code | Description | Unit | Price |
|--------------|--|--------|---------|
| ANS93-833-S | Gloves nitrile thin microflex blue small | 250pcs | \$21.00 |
| ANS93-833-M | Gloves nitrile thin microflex blue medium | 250pcs | \$21.00 |
| ANS93-833-L | Gloves nitrile thin microflex blue large | 250pcs | \$21.00 |
| ANS93-833-XL | Gloves nitrile thin microflex blue x-large | 250pcs | \$21.00 |



KONICA MINOLTA

SPAD 502 Plus chlorophyll meter

A lightweight and water resistant hand-held meter for lab and field, measuring the chlorophyll content of leaves without causing damage to the plant.

Measurements are taken by just closing the measuring head on the leaf. Since the leaf is not cut off or damaged in any other way, the same leaf can be measured as the plant grows.



| Item Code | Description | Unit | Price |
|-------------|---------------------------------|------|------------|
| MINA1RT-206 | SPAD-502 Plus chlorophyll meter | EA | \$3,180.00 |

Spectrophotometer CM-5

Stand alone full-featured spectrophotometer with maximum flexibility and versatility that covers all your applications

15% OFF LIST PRICE

Features and Benefits:

- Extremely easy operation with automatic calibration
- Top-port concept for Reflectance measurements as well as large transmission chamber for liquids and solids
- “Operation Wizard” with real time on screen user guidance and storage of user settings and data on USB stick



Compact Chroma meters / colorimeters

Konica Minolta’s color difference readers are portable, battery powered, and hand-held instrument for quick, precise color control.

| Item Code | Description | Unit | Price |
|-------------|--------------------------|------|------------|
| MINA80F-107 | Colour reader CR-10 Plus | EA | \$4,750.00 |
| MINA80F-217 | Colour reader CR-20 | EA | \$5,470.00 |

*Software included for CR-10 Plus and CR-20 readers



Stock up on lab bottles

Fisherbrand™ graduated laboratory bottles

- Polypropylene screw caps and pouring ring
- Borosilicate glass
- Autoclavable
- GL 45 screw thread
- Permanent white enamel marking spot



| Item Code | Description | Unit | Price |
|------------|--------------------------|------|----------|
| FSBFB33144 | Laboratory bottle 100 mL | EA | \$3.75 |
| FSBFB33145 | Laboratory bottle 250 mL | EA | \$3.95 |
| FSBFB33146 | Laboratory bottle 500 mL | EA | \$5.45 |
| FSBFB33147 | Laboratory bottle 1 L | EA | \$7.75 |
| FSBFB33148 | Laboratory bottle 2 L | EA | \$15.75 |
| FSBFB33180 | Laboratory bottle 5 L | EA | \$75.00 |
| FSBFB33181 | Laboratory bottle 10 L | EA | \$225.00 |
| FSBFB33182 | Laboratory bottle 20 L | EA | \$475.00 |

Time to clean up

Clean up, clear out, declutter and rejuvenate your bench and lab space. Let us help you out with our best cleaners for the job!



| Item Code | Description | Application | Unit | Price |
|-----------------|---|------------------------|------------------|----------|
| DEC90-5L | Decon 90 liquid concentrate | Manual cleaning | 5 L | \$85.00 |
| DECNEUTRACON-5L | Neutracon liquid | Manual cleaning | 5 L | \$85.00 |
| RBSRBS35 | RBS35 concentrated cleaning agent | | 5 L | \$100.00 |
| BSPSTE70/R | Ethanol denatured alcohol 70% sterile spray | General disinfection | 6 x 500 mL | \$70.00 |
| DIV12291 | Pyroneg liquid cleaner | Manual | 2 L | \$45.00 |
| DIV13231 | Pyroneg powder | Automatic cleaning | 3 kg | \$100.00 |
| LEI203330 | Virkon S | Virucidal disinfection | 50 x 5 g tablets | \$90.00 |
| LEI6569B | Virkon S broad spectrum | Virucidal disinfection | 1 kg | \$130.00 |

High quality Great value

Fisher Scientific™ accumet™ pH meters



Fisher Scientific accumet AE-150 pH Meter with Stand

Budget friendly, accurate and reliable basic pH meter

- Available as meter and stand only or in 3 different kits that include electrodes and temperature probes
- BNC input for electrode compatibility
- Automatic or manual temperature compensation
- pH slope and offset display
- One to five point pH calibration and mV offset adjustment for ORP
- 3 Years Warranty

| Item Code | Description | Unit | Price |
|----------------|--------------------------------|------|----------|
| ACC13636-AE150 | AE150 pH meter with stand | EA | \$539.00 |
| ACC13636-AE151 | AE150 pH meter - education set | EA | \$617.40 |
| ACC13636-AE152 | AE150 pH meter - bio set | EA | \$622.20 |
| ACC13636-AE153 | AE150 pH meter - 3 in one set | EA | \$632.30 |



Fisher Scientific accumet AB series meters

All models have the below features and are available as meter only or in kits that include electrodes and temperature probes

- AB-150: Able to measure pH/mV/Temperature
- AB-200: Able to measure pH/mV/Temperature/Conductivity/TDS/Resistivity/Salinity
- AB-250: Able to measure pH/mV/Temperature/Ion Selective Measurements
- 0.001 pH resolution and ± 0.002 pH accuracy
- Date and time logging to meet GLP requirements
- BNC input for electrode compatibility and RS-232 output for "21 CFR Part 11 Compliance"
- Password protection and 500 point memory with manual or automatic interval logging
- 3 Years Warranty

| Item Code | Description | Unit | Price |
|-----------------|-----------------------------------|------|------------|
| ACC13636-AB150A | accumet AB150 pH meter only | EA | \$830.60 |
| ACC13636-AB150 | accumet AB150 pH kit | EA | \$920.40 |
| ACC13636-AB150B | accumet AB150 pH bio kit | EA | \$984.50 |
| ACC13636-AB200A | AB200 pH/conductivity meter only | EA | \$1,252.50 |
| ACC13636-AB200 | accumet AB200 pH/conductivity kit | EA | \$1,440.80 |
| ACC13636-AB250A | accumet AB250 pH/ISE meter | EA | \$957.40 |
| ACC13636-AB250 | accumet AB250 pH/ISE kit | EA | \$1,010.80 |



Don't forget

We have a full range of chemicals,
check it out online at
www.thermofisher.com.au/chemicals

Don't forget

to check out more deals in our
lab essential promotion
www.thermofisher.com.au/labessentials

In Australia:

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, ordersau@thermofisher.com

Find out more at:

thermofisher.com.au/promotions

ThermoFisher
SCIENTIFIC

Offers valid till 31st March 2018. All offers subject to Thermo Fisher standard terms and conditions. Errors and omission exempted.
Specifications, terms and pricing are subject to change. Please consult your local sales representative for details. © 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. 1519099635